# CLEARFIL™ AP-X

LIGHT- CURED COMPOSITE RESIN

## ENGLISH INSTRUCTIONS FOR USE

### I. INTRODUCTION

CLEARFIL AP-X is a light-cured, radiopaque restorative composite resin which provides accurate color-matching, high polishability and excellent physical properties, making it ideal for both anterior and posterior restorations. It is formulated with optimal viscosity assuring easy handling and placement.

### II. INDICATIONS

CLEARFIL AP-X is indicated for the following restorative applications:

- Class I, II, V restorations of posterior teeth
- Class III, IV, V restorations of anterior teeth
- Cervical cavities or defects involving root surfaces

### III. CONTRAINDICATION

Patients with a history of hypersensitivity to methacrylate monomers

### IV. INCOMPATIBILITY

Do not use eugenol-containing materials for pulp protection or temporary sealing since the eugenol could retard the curing process.

### V. PRECAUTIONS

#### 1. Safety precautions

- Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomers.
- If any hypersensitivity, such as a rash or dermatitis, occurs from contact with the product, discontinue the use of the product and consult a physician.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers.
- Use caution to prevent the product from coming in contact with the skin or getting into the eyes. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from possible splashing material.
- If the product comes in contact with human body tissues, take the following actions.  
<If the product gets in the eyes>  
Immediately rinse the eye with copious amounts of water and consult a physician.  
<If the product comes in contact with the skin>  
Immediately wipe the material off using a cotton pledget (or gauze) moistened with alcohol and then wash with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Avoid looking directly at the curing light when curing the product.
- To prevent cross infection, avoid sharing the same paste among different patients. Dispense the paste on a mixing paper before use. Do not try to apply the paste directly from the syringe onto the patient's tooth.

#### 2. Handling and manipulation precautions

- If necessary, use a rubberdam in order to prevent moisture contamination.
- Do not use the product in conjunction with the other composite resin. Mixing materials may cause a change in physical properties, possibly a decrease, from the properties expected.
- If the adherent surface extends to uncut enamel, apply a phosphoric acid etching agent (e.g. K-ETCHANT GEL) and leave it in place for 10 seconds before washing and drying.
- Amalgam or other lining material remaining in the cavity will prevent the passage of light and the polymerization of the product. Completely remove any lining material when preparing the cavity.
- When light curing the product, note the depth of cure in this Instructions for Use.
- The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each

section separately.

- Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing unit guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light-evaluating device at appropriate intervals.
- The paste should be used as soon as practicable after being dispensed from the syringe. If the paste is to be left for a while before use, it should be covered with a light blocking plate.
- The product should be returned to room temperature before dispensing, if it has been taken out of a refrigerator. Failure to do this might cause breakage of the syringe.
- The use of the product is restricted to a licensed dental professional.

#### 3. Storage precautions

- Do not use after the expiration date. Note expiration date on the outside of package.
- Keep away from extreme heat or direct sunlight.
- The product must be stored at 2 - 25 °C / 36 - 77 °F when not in use.
- After the paste has been dispensed, the syringe should be capped securely as soon as possible to prevent the paste inside the syringe from being hardened by ambient light, and to keep foreign matter from entering the syringe.
- The product must be stored in proper places where only dental practitioners can access it.

### VI. SHADES SYSTEM AND COMPONENTS

#### 1. Shades

CLEARFIL AP-X is available in 14 shades separately; 11 shades correspond on a best-match basis with VITA Shade Guide plus 3 shades for special needs. Select an appropriate shade by matching the shade guide to the restoration site. The XL and HO shades offer a strong, lighter white color. The HO shade is translucent white and is suitable for sites requiring extra white color. The CL shade corresponds to A2 opaque and is suitable for cervical sites where opacity is required.

- Standard shades : A2, A3, A3.5, B2, B3, C1, C2, C3
- Cervical shades : A4, B4, C4 CL
- Translucent shades : XL, HO

#### 2. Components

Please see the outside of the package for contents and quantity.

#### 3. Ingredients

Principle ingredients;

- Bisphenol A diglycidylmethacrylate (Bis-GMA)
- Triéthylène glycol diméthacrylate (TEGDMA)
- Silanated barium glass filler
- Silanated silica filler
- Silanated colloidal silica
- dl-Camphorquinone

The total amount of inorganic filler is approx. 71 vol%. The particle size of inorganic fillers ranges from 0.02 to 17 µm.

### VII. CLINICAL PROCEDURES

- Cleaning tooth structure  
Be sure the cavity is adequately cleaned. An adequately cleaned cavity assures maximum adhesive performance.
- Moisture control  
Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.
- Cavity preparations  
Remove any infected dentin and prepare the cavity in the usual manner.
- Shade selection  
Select an appropriate shade using the VITA Lumin-VACUUM shade guide.

#### 5. Pulp protection

Any actual or near pulp exposure should be covered with a hard setting calcium hydroxide material. However, usually a cement lining or basing is not necessary. Do not use eugenol materials for pulp protection.

#### 6. Acid etching uncut enamel

If the resin restorative material will extend to uncut enamel, apply etching agent (e.g. K-ETCHANT GEL) to the enamel, let it stay for 10 seconds, wash with water, and then dry.

#### [CAUTION]

The use of a bonding agent alone does not adequately condition uncut enamel. Overfilling of resin onto unetched, uncut enamel could cause marginal discoloration.

#### 7. Application of bonding agent

Tooth surface treatment and bonding should be performed according to Instructions for Use of the bonding system used (e.g. CLEARFIL SE BOND or CLEARFIL TRI-S BOND).

#### 8. Placement and light curing of CLEARFIL AP-X

##### 1) Dispensing

Dispense the necessary amount of CLEARFIL AP-X. After dispensing the resin, turn the plunger counter-clockwise a half turn to prevent excess resin from escaping. Replace the syringe cap immediately to prevent premature setting of the resin. Use the resin within 5 minutes of dispensing. Extended exposure to the operating light will allow the resin to cure.

##### 2) Placement

Incremental placement and light-curing each increment is strongly recommended especially in deep cavities and in Class II cavities.

##### 3) Curing

Light-cure the resin with a dental visible light curing unit\* following the table. Hold the light tip as close to the resin as possible.

Table: Relation between curing time and depth of cure for each dental curing unit.

Type	Curing time (sec)	Depth of cure (mm)	
		XL, A2, A3, A3.5, B2, B3, C1, C2	A4, B4, C3, C4, CL, HO
Conventional halogen*	20	2.0	1.5
	40	2.0	2.0
Fast halogen*	5	2.0	1.5
	10	2.0	2.0
Plasma arc*	5	2.0	1.5
	10	2.5	2.0
LED*	20	2.0	1.5
	40	2.0	2.0

\*Dental curing unit

Type	Light source	Wavelength range and light intensity
Conventional halogen	Halogen lamp	Light intensity <sup>2)</sup> of 150 – 550 mW/cm <sup>2</sup> in wavelength range from 400 – 515 nm
Fast halogen	Halogen lamp	Light intensity <sup>2)</sup> of more than 550 mW/cm <sup>2</sup> in wavelength range from 400 – 515 nm
Plasma arc	Xenon lamp	Light intensity <sup>3)</sup> of more than 2000 mW/cm <sup>2</sup> in wavelength range from 400 – 515 nm, and light intensity of more than 450 mW/cm <sup>2</sup> in wavelength range from 400 – 430 nm
LED	Blue LED <sup>1)</sup>	Light intensity <sup>2)</sup> of more than 300 mW/cm <sup>2</sup> in wavelength range from 400 – 515 nm

1) Peak of emission spectrum: 450 – 480 nm

2) Evaluated according to ISO 10650-1.

3) Wavelength distribution and light intensity values measured with a spectro-radiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp

#### 9. Finishing

Contour the restoration and adjust the occlusion using a fine diamond point. Polish with silicon rubber points or polishing discs.

#### [CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

#### [WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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CLEARFIL, CLEARFIL AP-X, CLEARFIL TRI-S BOND and SE BOND are trademarks of KURARAY CO., LTD. VITA is a trademark of VITA Zahnfabrik, H. Rauter GmbH & Co. KG.
Fabricado por **Kuraray Noritake Dental Inc.** 1621 Sakazu, Kurashiki, Okayama 710-0801, Japan