

DENTAL LIGHT-CURE CORE BUILD-UP MATERIAL

CLEARFIL™ PHOTO CORE PLT

Flow Chart Sheet

Case A Core build-up of vital tooth

Case B Core build-up of non-vital tooth

- 1) Cementing the post using PANAVIA™ F 2.0
- 2) Cementing the post using another resin luting cement

S365EU-01



KURARAY MEDICAL INC.
1621 Sakazu, Kurashiki,
Okayama 710-8622, Japan

11/2003

DENTAL LIGHT-CURE CORE BUILD-UP MATERIAL

CLEARFIL™ PHOTO CORE PLT

Case A Core build-up of vital tooth

A

<p>Apply PRIMER gently</p>	<p>Dry with mild air flow</p>	<p>Apply BOND</p>	<p>Air flow gently</p>
<p>Light-cure</p>	<p>Placement</p>	<p>Light-cure * Refer to the instructions for use about light curing time.</p>	<p>Prepare of the abutment tooth</p>

DENTAL LIGHT-CURE CORE BUILD-UP MATERIAL

CLEARFIL™ PHOTO CORE PLT

Case B Core build-up of non-vital tooth (1/2)

1) Cementing the post using PANA VIA™ F 2.0

When using PANA VIA™ F 2.0, it is not necessary to use CLEARFIL™ SE BOND.

B

<p>Cementing the post using PANA VIA™ F 2.0</p>	<p>Place build-up directly on PANA VIA™ F 2.0 light or self-cure</p>	<p>Light-cure * Refer to the instructions for use about light curing time.</p>	<p>Prepare of the abutment tooth</p>

DENTAL LIGHT-CURE CORE BUILD-UP MATERIAL

CLEARFIL™ PHOTO CORE PLT

Case B Core build-up of non-vital tooth (2/2)

2) Cementing the post using another resin luting cement

<p>Apply PRIMER gently to the post and tooth surface</p>	<p>Dry with mild air flow</p>	<p>Apply BOND</p>	<p>Air flow gently</p>
<p>Light-cure</p>	<p>Build-up core</p>	<p>Light-cure * Refer to the instructions for use about light curing time.</p>	<p>Prepare of the abutment tooth</p>

B