

001 1562R402R




SELF-ADHESIVE RESIN CEMENT

PANAVIA™ SA Cement Plus Handmix
Flow Chart SheetKuraray Noritake Dental Inc.
1621 Sakazu, Kurashiki,
Okayama 710-0801, Japan**Indications for Use**

- [1] Cementation of crowns, bridges, inlays and onlays
- [2] Cementation of prosthetic restorations on implant abutments and frames*
- [3] Cementation of adhesion bridges and splints*
- [4] Cementation of posts and cores
- [5] Amalgam bonding*



* Refer to the IFU for [2], [3] and [5] of Indications.

**Indication 1 Cementation of Crowns**

Clean and dry the tooth surface, and then trial fit the prosthetic restoration.

1 Conditioning the prosthetic restoration

Follow the Instructions for Use of the restoration material. In the absence of specific instructions, we recommend the following procedure:

**2 Dispense an equal amount of Paste A & B. Mix Paste A & B for 10 seconds.***

* As necessary, blast with alumina powder, then ultrasonic clean and dry. The air pressure should be properly adjusted to suit the material and/or shape of the prosthetic restoration, using caution to prevent chipping.

* Refer to table 1 for working time.

3 Apply the mixed paste to the prosthetic restoration.*


* Refer to table 1 for working time.

4 Place the crown.**5 Light-cure for 2 to 5 seconds or chemical-cure for 2 to 4 minutes, then remove the excess cement.****6 Maintain isolation for 5 minutes.***

* For a translucent restoration, light-cure. Refer to table 2.

Indication 4 Cementation of Posts


Clean and dry the cavity, and then trial fit the post.

1 Dispense an equal amount of Paste A & B. Mix Paste A & B for 10 seconds.*

* Refer to table 1 for working time.

2 Apply over the entire adherend surface of the post, or the entire tooth surface within the cavity.*

* Refer to table 1 for working time.

3 Place the post quickly into the cavity, slightly vibrating it to prevent air bubbles from entering the root canals.**4 Spread the excess paste over the coronal base and post head. Light-cure the margins of the post.***

* Refer to table 2.

5 Place the core buildup composite resin.*

* Prepare the abutment tooth 10 minutes after placing the dental post.

Table 1 : Working time

Working time after mixing (23°C/ 73°F)	2 min.
Working time after insertion of the paste into the cavity (37°C/ 99°F)	40 sec.

Table 2 : Curing time for type of light source

Type	Light source	Light Intensity	Curing time
High-intensity BLUE LED	BLUE LED	More than 1500 mW/cm ²	Twice for 3 to 5 sec.
BLUE LED		800-1400 mW/cm ²	10 sec.
Halogen	Halogen lamp	More than 400 mW/cm ²	10 sec.

* For the light intensity, refer to the IFU of the dental curing unit.