DENTAL DUAL-CURED ADHESIVE RESIN CEMENT

PANAVIA™ F 2.0 Complete Kit

[ENGLISH] INSTRUCTIONS FOR USE

I. INTRODUCTION

PANAVIA F 2.0 Complete Kit consists of PANAVIA F 2.0 Paste, PANAVIA F 2.0 ED PRIMER II, CLEARFIL CERAMIC PRIMER, ALLOY PRIMER and PANAVIA F 2.0 OXYGUARD II. PANAVIA F 2.0 Paste is a dual-cure (light- and/or self-cure), radiopaque resin-based cement for conventional porcelain, ceramic, hybrid ceramics (e.g. ESTENIA C&B), composite resin, and metal restorations. PANAVIA F 2.0 ED PRIMER II is a tooth surface conditioner and consists of Liquid A and Liquid B. CLEARFIL CERAMIC PRIMER is a silane-coupling agent that provides an enhanced adhesive surface to conventional porcelain, ceramic, hybrid ceramics or composite resin. ALLOY PRIMER improves the bond strength for precious alloy and PANAVIA F 2.0 Paste. PANAVIA F 2.0 OXYGUARD II is an oxygen-blocking agent to allow the PANAVIA F 2.0 Paste to polymerize when not light-cured.

II. INDICATIONS

PANAVIA F 2.0 Complete Kit is indicated for the following uses:

- Cementation of crowns, bridges, inlays and onlays made of conventional porcelain, ceramic, hybrid ceramics, composite resin or metal
- Cementation of veneers
- Cementation of adhesion bridges
- Cementation of metal cores, resin cores, metal posts or glass-fiber posts
- Amalgam bonding

[NOTE]

Use a cement shade appropriate to the individual case.

The shades of adhesive cements and applicable cases:

Restoration	Cement shade	TC, Light	White	Opaque
Crown, bridge, inlay or onlay made of metal		○	○	○
Crown, bridge, inlay or onlay made of conventional porcelain, ceramic, hybrid ceramics or composite resin		○	△	△
Veneer		○	△	△
Core or post		○	○	○
Adhesion bridge and splint	anterior	△	○	○
	posterior	○	○	○
Bonded amalgam restoration		○	○	○

○: RECOMMENDED ○: ADEQUATE △: NOT RECOMMENDED

III. CONTRAINDICATION

- Patients with a history of hypersensitivity to methacrylate monomers
- Patients with a history of hypersensitivity to acetone
- Patients with a history of hypersensitivity to acid

IV. POSSIBLE SIDE EFFECTS

The oral mucosal membrane may turn a whitish color when contacted by PANAVIA F 2.0 ED PRIMER II, CLEARFIL CERAMIC PRIMER or ALLOY PRIMER due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. In individual cases ulceration has been reported.

V. INCOMPATIBILITIES

- Do not use eugenol-based materials for pulp protection or temporary sealing, since the eugenol could retard the curing process.
- Do not use hemostatic agents containing ferric compounds, since these materials may impair adhesion and may cause discoloration at the tooth margin or surrounding gingiva, due to ferric ions which may remain.

VI. PRECAUTIONS

- Safety precautions**
 - This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
 - If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
 - Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
 - If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye> Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin or the oral mucosa> Immediately wipe the area with a cotton pellet or a gauze moistened with alcohol, and rinse with copious amounts of water.
 - Use caution to prevent the patient from accidentally swallowing the product.
 - Do not use the same disposable brush tip for different patients to prevent cross contamination. Discard the tip after use and sterilize the brush tip handle after each patient.
 - Avoid direct contact with the skin and / or soft tissue to prevent hypersensitivity. Wear gloves or take appropriate precautions when using the product.
 - If a patient or operator feels sick from inhaling the acetone contained in the product, allow them to rest and breathe fresh air.
 - When using with preformed stainless posts, the post should not contact metal restorations. Cover the post with composite resin.

2. Handling and manipulation precautions

- [PANAVIA F 2.0 Paste]
 - The paste should be brought to room temperature for 15 minutes or more after it is removed from the refrigerator; this will restore the normal viscosity of the paste. In addition, this helps to prevent water contamination from possible refrigerator moisture.
 - The paste must be used within 3 minutes after mixing.
 - The paste should not be applied to a restoration surface that has been treated with PANAVIA F 2.0 ED PRIMER II; otherwise the paste will cure on contact with the PANAVIA F 2.0 ED PRIMER II. The procedure should be carried out within the working time listed in the table, "Working time of PANAVIA F 2.0 Paste" in VIII. CLINICAL PROCEDURES.
 - The paste contains a light-cure catalyst that is highly photo-reactive. Use the light blocking plate to avoid exposing the material to an operating light or natural light (sunlight from windows). During curation, adjust the angle and/or distance of the dental operating light to reduce the intensity of light entering the oral cavity and prevent further polymerization of the paste.
 - Do not use a lentulo spiral for loading the PANAVIA F 2.0 Paste into the root canal. If the paste is loaded into the root canal using a composite resin syringe, the polymerization of the paste is accelerated. It is necessary to fit the post as quickly as possible.
 - If you want to place dental posts into multiple root canals of a posterior tooth, complete the post placement of one root canal before proceeding with another, and prevent the excess paste from entering the other root canals.
 - Any paste that has hardened on the gingiva will come off; however, try to remove the excess paste before it sets whenever possible.
 - The Opaque paste should not be light-cured, but allowed to cure by using PANAVIA F 2.0 OXYGUARD II. It has a low curing depth.
- [PANAVIA F 2.0 ED PRIMER II]
 - ED PRIMER II should be brought to room temperature before using. After taking out of the refrigerator, the product must be left standing for 15 minutes, until it comes to room temperature; otherwise, bubbles will form in the liquid when dispensing or it will ozone after use.
 - Dispense ED PRIMER II Liquid A and B from each container by pointing the container downward as vertically as possible; careful dispensing is necessary to assure equal amounts of both liquids otherwise, the physical property of the product may be affected.
 - ED PRIMER II Liquid A and B should be mixed whenever used. Do not use them separately or individually as a single agent.
 - Use the mixture of ED PRIMER II Liquid A and B as soon as possible after mixing. The mixture must be used within 5 minutes after mixing.
 - When luting to a precious metal abutment, apply ALLOY PRIMER before applying ED PRIMER II. When luting to an abutment made of a resin containing inorganic filler, apply CLEARFIL CERAMIC PRIMER before applying ED PRIMER II.
 - If the treated surface is contaminated with saliva or tissue exudates, wash it with water, dry, or clean with ethanol, and treat with ED PRIMER II again; otherwise, optimal adhesion will be impaired.

[CLEARFIL CERAMIC PRIMER]

- CLEARFIL CERAMIC PRIMER should be brought to room temperature before using. After taking out of the refrigerator, the product must be left standing for 15 minutes, until it comes to room temperature; otherwise, bubbles will form in the liquid when dispensing or it will ozone after use.
- Do not mix CLEARFIL CERAMIC PRIMER with bonding agents or use it in conjunction with other product; intermixing materials may cause a change in physical properties, including a possible decrease from the properties expected.

4. The necessary amount of paste for a typical application is:

Number of rotations of syringe	Applications
Half a turn	Inlays and onlays
1 turn	Crowns

[NOTE]

- If the paste is dispensed by turning the plunger a quarter of a turn, the performance of the product could be impaired when the paste hardens.
- If not used immediately, it should be covered with a light blocking plate.

2) Mixing A Paste and B Paste

- Use sufficient A Paste and B Paste on the paper pad for 20 seconds. Be sure there is no water mist on the paper pad or spatula before using them, the presence of water mist could shorten the working time of the mixed paste.
- Perform cementation soon after treating the restoration surface with CLEARFIL CERAMIC PRIMER.

- Do not use CLEARFIL CERAMIC PRIMER near an open flame. It contains ethanol, a flammable substance.
- If the treated surface is contaminated with saliva or tissue exudates after CLEARFIL CERAMIC PRIMER is applied, clean with a phosphoric acid (e.g. K-ETCHANT GEL), wash with water, dry and apply CLEARFIL CERAMIC PRIMER again.

Working time of PANAVIA F 2.0 Paste	Working steps	Working time
1.	After dispensing the pastes (by turning the syringe the same amount of turn)	15 min.
2.	After mixing the pastes (for 20 sec.)	3 min.
3.	Placing the restoration under pressure --In the case of root canal	60 sec. 40 sec.

[ALLOY PRIMER]

- If the treated surface is contaminated with saliva or tissue exudates after ALLOY PRIMER is applied, clean the restoration surface with a cotton pellet moistened with ethanol and apply the product again.
- Do not use ALLOY PRIMER near an open flame. It contains acetone, a flammable substance.

[Dental light curing unit]

- Do not look directly at the light source. Protective glasses are recommended.
- If the light curing unit has low intensity, poor curing will result. Check the service life of the lamp and examine the emitting tip for contamination at periodic intervals. It is also advisable to check the time required to cure PANAVIA F 2.0 Paste by light-curing a sample of PANAVIA F 2.0 Paste mix prior to starting the treatment (36 - 46°F).
- Check the conditions required to cure the paste mix by referring to the light-curing times listed in this Instructions for Use before using the product.

The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

[Common precautions]

- Provide moisture and contamination control with the use of a rubber dam.
- Any exposed pulp or areas near the pulp should be covered with a hard setting calcium hydroxide material.
- Clean the cavity sufficiently to prevent poor bonding. If the restoration surface is contaminated with saliva or tissue exudates, wash it thoroughly and dry before cementation.
- To prevent poor performance and handling characteristics, observe the specified light-curing times and other handling requirements.
- Do not mix the product with any other dental material.
- Be careful not to cut your fingers on the sharp edges of the instruments.
- Do not use the product for any purposes other than those specified in restorations.
- The use of this product is restricted to a licensed dental professional.

3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- PANAVIA F 2.0 Complete Kit must be stored in a refrigerator (2 - 8°C / 36 - 46°F) when not in use. All the components except for ALLOY PRIMER must be stored in a refrigerator (2 - 8°C / 36 - 46°F) when not in use. ALLOY PRIMER must be stored at 2-25°C/36-77°F when not in use.
- Keep away from extreme heat or direct sunlight.
- The bottle or syringe cap should be replaced as soon as the liquid or resin has been dispensed from the bottle or syringe. This prevents evaporation of volatile ingredients.
- The product must be stored in a proper place where only dental practitioners can access.

VII. COMPONENTS

1. Shades

PANAVIA F 2.0 Paste is available in 4 shades;

TC, Light, White or Opaque

2. Components

Please see the outside of the package for contents and quantity.

- PANAVIA F 2.0 Paste (A Paste / B Paste) : TC, Light, White or Opaque
- PANAVIA F 2.0 ED PRIMER II (Liquid A/Liquid B)
- CLEARFIL CERAMIC PRIMER
- ALLOY PRIMER
- PANAVIA F 2.0 OXYGUARD II
- Accessories
 - Spatula
 - Mixing pad
 - Disposable brush tips
 - Brush tip handle
 - Mixing dish
 - Light blocking plate
 - Disposable nozzles

3. Ingredients

1) PANAVIA F 2.0 Paste (A Paste/B Paste)

Principal ingredients

- A Paste
 - 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
 - Hydrophobic aromatic dimethacrylate
 - Hydrophobic aliphatic dimethacrylate
 - Hydrophilic aliphatic dimethacrylate
 - Silanated silica filler
 - Silanated colloidal silica
 - Camphorquinone
 - Catalysts
 - Initiators
- B Paste
 - Hydrophobic aromatic dimethacrylate
 - Hydrophobic aliphatic dimethacrylate
 - Hydrophilic aliphatic dimethacrylate
 - Silanated barium glass filler
 - Surface treated sodium fluoride
 - Catalysts
 - Accelerators
 - Pigments

The total amount of inorganic filler is approx. 59 vol%. The particle size of inorganic fillers ranges from 0.04 μm to 19 μm.

2) PANAVIA F 2.0 ED PRIMER II

Principal ingredients

- Liquid A
 - 2-Hydroxyethyl methacrylate (HEMA)
 - 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
 - Water
 - N-Methacryloyl-5-aminosilicic acid (5-NMSA)
 - Accelerators
- Liquid B
 - N-Methacryloyl-5-aminosilicic acid (5-NMSA)
 - Water
 - Catalysts
 - Accelerators

The total amount of inorganic filler is approx. 59 vol%. The particle size of inorganic fillers ranges from 0.04 μm to 19 μm.

3) CLEARFIL CERAMIC PRIMER

Principal ingredients

- Ethanol
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- 3-Methacryloyloxypropyl trimethoxysilane

4) ALLOY PRIMER

Principal ingredients

- Acetone
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- 6-(4-Vinylbenzyl-N-propylamino-1,3,5-triazine-2,4-dithione)
 - Principal ingredients
 - Glycerol
 - Polyethylene glycol
 - Catalysts
 - Accelerators
 - Dyes

VIII. CLINICAL PROCEDURES

Preparing PANAVIA F 2.0 Paste

1) Dispensing the pastes

- Align markings on the nut with the reference line on the plunger and turn the syringe to dispense the necessary amount of paste. A minimum turn of the syringe should be half a turn.
- Equal amounts of A Paste and B Paste should be dispensed.
- The amount of paste dispensed from a 45° rotation of the syringe could be inaccurate; therefore, discard the syringe before using the last portion.

Table: Dental curing unit

Type	Light source	Wavelength range and light intensity
Conventional halogen lamp	Halogen lamp	Light intensity ¹⁾ of 300 - 550 mW/cm ² in wavelength range from 400 - 515 nm
Fast halogen lamp	Halogen lamp	Light intensity ¹⁾ of more than 550 mW/cm ² in wavelength range from 400 - 515 nm
Plasma arc	Xenon lamp	Light intensity ¹⁾ of more than 2000 mW/cm ² in wavelength range from 400 - 515 nm, and light intensity of more than 450 mW/cm ² in wavelength range from 400 - 430 nm
LED	Blue LED ¹⁾	Light intensity ²⁾ of more than 300 mW/cm ² in wavelength range from 400 - 515 nm

- Peak of emission spectrum: 450 - 480 nm
- Evaluated according to ISO 10665-1.
- Wavelength distribution and light intensity values measured with a spectro-radiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp

Table: Light curing time

Dental curing unit	Light curing time
Conventional halogen LED	20 seconds
Fast halogen lamp	5 seconds

[CAUTION]

- The Opaque paste should not be light-cured; allow to cure by using PANAVIA F 2.0 OXYGUARD II. It has a low curing depth.
- Self-curing method:
 - Use PANAVIA F 2.0 OXYGUARD II to cure the mixed paste as follows: With a disposable brush tip, apply OXYGUARD II to the margin. After 3 minutes, remove OXYGUARD II with a cotton roll and water spray.

A-10. Finishing

Remove the excess paste adhered to both surface by polishing.

B. Standard procedure III [indications 4]: for core build up

[4] Cementation of metal cores, resin cores, metal posts or glass-fiber posts

[NOTE]

This procedure is for use with a preformed post and composite resin core build-up. For the cementation of cores, refer to standard procedure I (Working time of PANAVIA F 2.0 Paste is 40 seconds in the case of root canal. See table "Working time of PANAVIA F 2.0 Paste").

B-1. Cavity preparation and trial fit of the post

- Prepare the endodontically filled root canals for post placement in the usual manner. Provide moisture and contamination control with a rubber dam.
- Trial fit a post of appropriate length and thickness in the prepared cavity. Cut and trim the post as necessary. Wipe away any contamination from the surface of the post using a piece of gauze or a cotton pad soaked with ethanol.

B-2. Post surface treatment

For metal post

- If the post you want to use is not yet surface treated, sandblast it. Refer to the instruction for use of the specific post.
- Clean the post in an ultrasonic unit for 2 minutes.
- Apply ALLOY PRIMER to the cleaned post with a disposable brush tip.

[NOTE]

There is no need to treat the adherent surface if it is non-precious metal.

For glass-fiber post

- Apply phosphoric acid (e.g. K-ETCHANT GEL) to the post surface with a disposable brush tip, leave it for 5 seconds, then wash with water and dry the surface.
- Apply CLEARFIL CERAMIC PRIMER to the post surface with a disposable brush tip.

[CAUTION]

Do not sandblast glass-fiber posts; this could damage the post fibers. Damage to the fibers degrades the strength to adhesive resin materials. CLEARFIL CERAMIC PRIMER contains volatile ethanol. Dry the post surface sufficiently by using high-pressure, oil-free air flow; otherwise, the adhesion strength will be impaired. Observe the described drying method and treatment time to ensure optimum adhesion.

B-3. Application of PANAVIA F 2.0 ED PRIMER II to the root canal and the cavity wall

For the application of ED PRIMER II to the root canal and the cavity wall, refer to A-5.

B-4. Preparing PANAVIA F 2.0 Paste

Refer to "Preparing PANAVIA F 2.0 Paste" in "VIII CLINICAL PROCEDURES".

B-5. Placing the post

- Apply the mixed PANAVIA F 2.0 Paste to the post.
- After applying the mixed paste to the post, insert it into the root canal quickly. It is advisable to lightly vibrate the post while inserting it into the root canal to prevent the entrapping of air bubbles. Cementation should be completed within 40 seconds.

[CAUTION]

- Do not use a lentulo spiral for loading the PANAVIA F 2.0 Paste into the root canal. If the paste is loaded into the root canal using a composite resin syringe, the polymerization of the paste is accelerated. It is necessary to fit the post as quickly as possible.
- If multiple posts are to be fitted to a single tooth, use caution to prevent excess PANAVIA F 2.0 Paste from flowing into other root canals.

B-6. Spreading the excess PANAVIA F 2.0 Paste

Using a disposable brush tip, spread the excess paste over the coronal base and post head.

B-7. Light-curing

Light-cure PANAVIA F 2.0 Paste on the remaining crown and post head. See table "Dental curing unit" and "Light-curing time in A-9."

[CAUTION]

For the application of Opaque paste, light-cure is not necessary. Please proceed to the next step.

B-8. Core build-up

After cementing the post in the canal, place a core build-up composite resin (e.g. CLEARFIL PHOTO CORE or CLEARFIL DC CORE AUTOMIX) for preparation of the abutment tooth according to the manufacturer's instructions.

C. Standard procedure III [indications 5]

5) Amalgam bonding

C-1. Cleaning of tooth structure

Clean the cavity and provide moisture control in the usual manner.

C-2. Application of PANAVIA F 2.0 ED PRIMER II to the cavity

For the application of ED PRIMER II to the cavity, refer to A-5.

C-3. Preparing PANAVIA F 2.0 Paste

Refer to "Preparing PANAVIA F 2.0 Paste" in "VIII CLINICAL PROCEDURES".

C-4. Placing the amalgam

- Apply the mixed PANAVIA F 2.0 Paste to the cavity.
 - Apply a thin, even layer of the mixed paste to the entire cavity surface primed with PANAVIA F 2.0 ED PRIMER II being careful to avoid entrapping air.

[CAUTION]

Because PANAVIA F 2.0 ED PRIMER II accelerates the set of PANAVIA F 2.0 Paste, the paste should be applied to the primed cavity quickly.

C-5. Amalgam filling

The triturated amalgam should be condensed on the unset PANAVIA F 2.0 Paste. Occlusal carving can be accomplished in the normal manner.

C-5. Removing the excess PANAVIA F 2.0 Paste

Any excess paste remaining at the margins can be removed with a disposable brush tip.

C-6. Curing PANAVIA F 2.0 Paste

For curing the mixed paste along the cement margin, refer to A-9.

C-7. Finishing

Remove the excess PANAVIA F 2.0 Paste adhered to tooth structure by polishing.

[CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

[WARNING]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, indirect or special, arising out of the application or use of, or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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CLEARFIL, PANAVIA, OXYGUARD, CLEARFIL PHOTO CORE and CLEARFIL DC CORE AUTOMIX are trademarks of KURARAY CO., LTD.