## DUAL CURE DENTAL ADHESIVE SYSTEM

# PANAVIA F 2.0

### ENGLISH INSTRUCTIONS FOR USE

#### I. INTRODUCTION

PANAVIA F 2.0 consists of PANAVIA F 2.0 Paste, ED PRIMER II and OXYGUARD II. PANAVIA F 2.0 Paste is a dual-cure light- and/or self-cure), radiopaque resin-based cement for porcelain, composite resin, and metal restorations. ED PRIMER II is a tooth surface conditioner and consists of Liquid A and Liquid B. OXYGUARD II is an oxygen-blocking agent to allow the PANAVIA F 2.0 Paste to polymerize when not light-cured.

#### II. INDICATIONS

PANAVIA F 2.0 is indicated for the following uses:

- Cementation of metal crowns and bridges, inlays and onlays
- Cementation of porcelain crowns, inlays, onlays and veneers
- Cementation of composite resin crowns, inlays, and onlays
- Cementation of adhesion bridges
- Cementation of endodontic cores and prefabricated posts
- Amalgam bonding

#### [NOTE]

Use a cement shade appropriate to the individual case.

The shades of adhesive cements and applicable cases:

Restoration	Cement shade	TC, Light	White	Opaque
Metal inlays and onlays; metal crowns and bridges		○	○	○
Porcelain or composite inlays, onlays and veneers		○	△	△
Preformed posts and cast metal cores		○	○	○
Adhesion bridges and splint	anterior	△	○	○
	posterior	○	○	○
Bonded amalgam restorations		○	○	○

○: RECOMMENDED ○: ADEQUATE △: NOT RECOMMENDED

#### III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers

#### IV. POSSIBLE SIDE EFFECTS

The oral mucosal membrane may turn a whitish color when contacted by ED PRIMER II due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. In individual cases, ulceration has been reported.

#### V. INCOMPATIBILITIES

- Do not use eugenol-based materials for pulp protection or temporary sealing, since the eugenol could retard the curing process.
- Do not use hemostatic agents containing ferric compounds, since these materials may impair adhesion and may cause discoloration at the tooth margin or surrounding gingiva, due to ferric ions which may remain.

#### VI. PRECAUTIONS

##### 1. Safety precautions

- This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
- If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
- If the product comes in contact with human body tissues, take the following actions:
  - <If the product gets in the eye>
    - Immediately wash the eye with copious amounts of water and consult a physician.
    - <If the product comes in contact with the skin or the oral mucosa>
      - Immediately wipe the area with a cotton pellet or a gauze moistened with alcohol, and rinse with copious amounts of water.
  - Use caution to prevent the patient from accidentally swallowing the product.
  - Do not use the same disposable brush tip for different patients to prevent cross contamination. Discard the tip after use and sterilize the brush tip handle after each patient.
  - Avoid direct contact with the skin and / or soft tissue to prevent hypersensitivity. Wear gloves or take appropriate precautions when using the product.
  - If a patient or operator feels sick from inhaling the organic solvents contained in the product, allow them to rest and breathe fresh air.
  - When using with preformed stainless posts, the post should not contact metal restorations. Cover the post with composite resin.

##### 2. Handling and manipulation precautions

[PANAVIA F 2.0 Paste]

- The paste should be brought to room temperature for 15 minutes or more after it is removed from the refrigerator; this will restore the normal viscosity of the paste. In addition, this helps to prevent water contamination from possible refrigerator moisture.
- The paste must be used within 3 minutes after mixing.
- The paste should not be applied to a restoration surface that has been treated with ED PRIMER II; otherwise the paste will cure on contact with ED PRIMER II. The procedure should be carried out within the working time listed on the table. "Working time of PANAVIA F 2.0 Paste" in VIII. CLINICAL PROCEDURES.
- The paste contains a light-cure catalyst that is highly photo-reactive. Use the light blocking plate to avoid exposing the material to an operating light or natural light (sunlight from windows). During cementation, adjust the angle and/or distance of the dental operating light to reduce the intensity of light entering the oral cavity and prevent premature polymerization of the paste.
- Do not use a lentulo spiral for loading the PANAVIA F 2.0 Paste into the root canal. If the paste is loaded into the root canal using a composite resin syringe, the polymerization of the paste is accelerated. It is necessary to fit the post as quickly as possible.
- If you want to place dental posts into multiple root canals of a posterior tooth, complete the post placement of one root canal before proceeding with another, and prevent the excess paste from entering the other root canals.
- Any paste that has hardened on the gingiva will come off; however, try to remove the excess paste before it sets, whenever possible.
- The Opaque paste should not be light-cured, but allowed to cure by using OXYGUARD II. It has a low curing depth.

#### [ED PRIMER II]

- ED PRIMER II should be brought to room temperature before using. After taking out of the refrigerator, the product must be left standing for 15 minutes, until it comes to room temperature; otherwise, bubbles will form in the liquid when dispensing or it will ooze after use.
- Dispense ED PRIMER II Liquid A and B from each container by pointing the container downward as vertically as possible; careful dispensing is necessary to assure equal amounts of both liquids otherwise, the physical property of the product may be affected.
- ED PRIMER II Liquid A and B should be mixed whenever used. Do not use them separately or individually as a single agent.
- Use the mixture of ED PRIMER II Liquid A and B as soon as possible after mixing. The mixture must be used within 5 minutes after mixing.
- When luting to a precious metal abutment, apply metal adhesive primer (e.g. ALLOY PRIMER) before applying ED PRIMER II. When luting to an abutment made of a resin containing inorganic filler, apply silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER or the mixture of CLEARFIL PORCELAIN BOND ACTIVATOR and CLEARFIL SE BOND PRIMER or CLEARFIL TRI-S BOND) before applying ED PRIMER II.
- If the treated surface is contaminated with saliva or tissue exudates, wash it with water, dry, or clean with ethanol, and treat with ED PRIMER II again; otherwise, optimal adhesion will be impaired.

#### [Dental light curing unit]

- Do not look directly at the light source. Protective glasses are recommended.
- If the light curing unit has low intensity, poor curing will result. Check the service life of the lamp and examine the emitting tip for contamination at periodic intervals. It is also advisable to check the time required to cure PANAVIA F 2.0 Paste by light-curing a sample of PANAVIA F 2.0 Paste mixture prior to starting the treatment.

#### [CAUTION]

- The working time of mixed paste could vary if the mixing is insufficient.
- The paste must be used within 3 minutes after mixing. When the mixed paste is left more than 3 minutes after mixing, the working time of PANAVIA F 2.0 Paste is shorter than that described on the above table. Use the mixed paste as soon as possible after dispensing and mixing.

#### A. Standard procedure I (Indications [1] to [4]; for cementation)

- Cementation of metal crowns and bridges, inlays and onlays
- Cementation of porcelain crowns, inlays, onlays and veneers
- Cementation of composite resin crowns, inlays, and onlays
- Cementation of adhesion bridges

#### A-1. Cleaning the cavity or abutment surface

- Remove the temporary sealing material and temporary cement in the usual manner; clean the cavity and provide moisture control and contamination control.
- Try-in the prosthetic restoration to check the fit in the cavity or abutment.

#### A-2. Sandblasting the prosthetic restoration surface

For metal, porcelain or composite resin restorations

As necessary, roughen the adhering surface by sandblasting with 50 to 50 μm alumina powder at an air pressure of 0.1 to 0.4 MPa (14-5.8 PSI). The air pressure should be properly adjusted to suit the material and/or shape of the prosthetic restoration, using caution to prevent chipping. After sandblasting, clean the prosthetic restoration by using an ultrasonic unit for 2 minutes, followed by drying it with an air stream.

#### A-3. Pretreatment of the prosthetic restoration surface

For porcelain or composite resin restorations

- Apply phosphoric acid (e.g. K-ETCHANT GEL) to the adherent surface, leave it for 5 seconds, then wash the surface with water and dry the surface.
- Dispense the necessary amounts of silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER or the mixture of CLEARFIL PORCELAIN BOND ACTIVATOR and CLEARFIL SE BOND PRIMER or CLEARFIL TRI-S BOND) into a well of the mixing dish immediately before application. Apply the agent (e.g. CLEARFIL CERAMIC PRIMER or the mixture of CLEARFIL PORCELAIN BOND ACTIVATOR and CLEARFIL SE BOND PRIMER or CLEARFIL TRI-S BOND) to the internal surface of the restoration with a disposable brush tip. After application, dry the entire restoration surface sufficiently using mild, oil-free air flow.
- Dispense the necessary amounts of silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER or the mixture of CLEARFIL PORCELAIN BOND ACTIVATOR and CLEARFIL SE BOND PRIMER or CLEARFIL TRI-S BOND) into a well of the mixing dish immediately before application. Apply the agent (e.g. CLEARFIL CERAMIC PRIMER or the mixture of CLEARFIL PORCELAIN BOND ACTIVATOR and CLEARFIL SE BOND PRIMER or CLEARFIL TRI-S BOND) to the internal surface of the restoration with a disposable brush tip.

#### [NOTE]

After treating the restoration surface, proceed with cementation quickly.

#### For precious metal restorations

Apply metal adhesive primer (e.g. ALLOY PRIMER) to the restoration surface with a disposable brush tip.

#### [NOTE]

- After treating the restoration surface, proceed with cementation quickly.
- There is no need to treat the adhering surface if it is non-precious metal.

#### A-4. Pretreatment of the cavity or abutment surface

Treatment of enamel

When cementing to uncut enamel or using with adhesive bridges and porcelain veneers, apply phosphoric acid (e.g. K-ETCHANT GEL) to the enamel surface and leave it for 10 seconds, then wash the surface with water and dry the surface.

#### Treatment of precious metal surface

When precious metal is used, apply metal adhesive primer (e.g. ALLOY PRIMER) to the metal surface with a disposable brush tip.

#### A-5. Application of ED PRIMER II to the cavity or abutment surface

- Dispense one drop each of ED PRIMER II Liquid A and B into a well of the mixing dish and mix immediately before application.

#### [CAUTION]

Use the mixture of ED PRIMER II Liquid A and B as soon as possible after mixing. The mixture must be used within 5 minutes after mixing.

- Apply the mixed ED PRIMER II to the cavity or abutment surface with a disposable brush tip. Leave it in place for 30 seconds. Use caution not to allow saliva or tissue exudates to contact the treated surfaces.
- After conditioning the cavity and abutment surface for 30 seconds, use a cotton pellet to carefully remove any excess ED PRIMER II from the cavity or abutment, especially from the corners of the cavity and shoulders at the margins of the abutment.
- After wiping away the excess ED PRIMER II, thoroughly dry the entire cavity and abutment surface sufficiently by blowing mild, oil-free air. Use a vacuum aspirator to prevent the ED PRIMER II from dispersing.

#### [CAUTION]

- ED PRIMER II should be applied on the entire surface of the tooth structure or abutment surface. Do not apply it to the restoration.
- Dry ED PRIMER II completely with gentle air flow. Remember that a pool of ED PRIMER II will cause quick polymerization of the adhesive cement. Do not rinse.

#### A-6. Preparing PANAVIA F 2.0 Paste

Refer to the foregoing section "Preparing PANAVIA F 2.0 Paste " in this chapter.

#### A-7. Cementing the prosthetic restoration

- Apply the mixed PANAVIA F 2.0 Paste to the restoration.

#### [CAUTION]

- DO NOT apply PANAVIA F 2.0 Paste to the tooth surface primed with ED PRIMER II as this will accelerate the set of PANAVIA F 2.0 Paste.
- Cement the restoration to the cavity or abutment tooth. Cementation should be completed within 60 seconds.

#### A-8. Removing excess paste

Any excess paste remaining at the margins can be removed with a disposable brush tip.

#### A-9. Curing PANAVIA F 2.0 Paste

Cure the mixed PANAVIA F 2.0 Paste along the cement margin, using either of the following two methods:

##### (1) Light-curing method:

Light-cure the paste along the cement margin with a dental curing unit (See the table "Dental curing unit") for the specified length of time shown on the table "Light curing time".

Type	Light source	Wavelength range and light intensity
Conventional halogen	Halogen lamp	Light intensity* of 300 - 550 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm
Fast halogen	Halogen lamp	Light intensity* of more than 550 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm
Plasma arc	Xenon lamp	Light intensity* of more than 2000 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm, and light intensity of more than 450 mW/cm <sup>2</sup> in wavelength range from 400 - 430 nm
LED	Blue LED <sup>1)</sup>	Light intensity* of more than 300 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm

- Peak of emission spectrum: 450 - 480 nm
- Evaluated according to ISO 10650-1.
- Wavelength distribution and light intensity values measured with a spectro-radiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp
- Evaluated according to ISO 10650-2.

#### Table: Light curing time

Dental curing unit	Light curing time
Conventional halogen	20 seconds
LED	5 seconds
Plasma arc	

#### [CAUTION]

The Opaque paste should not be light-cured; allow to cure by using OXYGUARD II. It has a low curing depth.

#### (2) Self-curing method:

Use OXYGUARD II to cure the mixed paste as follows: With a disposable brush tip, apply OXYGUARD II to the margin. After 3 minutes, remove OXYGUARD II with a cotton roll and water spray.

#### A-10. Finishing

Remove the excess paste adhered to tooth surface by polishing.

#### B. Standard procedure II (indications [5] : for core build up)

- Cementation of endodontic cores and prefabricated posts

#### [NOTE]

This procedure is for use with a preformed post and composite resin core build-up. For the cementation of metal cores, refer to standard procedure I (The working time of PANAVIA F 2.0 Paste is 40 seconds in the case of root canal. See the table "Working time of PANAVIA F 2.0 Paste" ).

#### B-1. Cavity preparation and trial fit of the metal post

- Prepare the endodontically filled root canals for post placement in the usual manner. Provide moisture control and contamination control with a rubber dam.

(2) Trial fit a metal post of appropriate length and thickness in the prepared cavity. Cut and trim the post as necessary. Wipe any wax or contamination from the surface of the post using a piece of gauze or a cotton pad soaked with ethanol.

#### B-2. Post surface treatment

If the post you want to use is not yet surface treated, sandblast it. Refer to the instruction for use of the specific post.

- Clean the post in an ultrasonic unit for 2 minutes.
- Apply the metal adhesive primer (e.g. ALLOY PRIMER) to the cleaned post with a disposable brush tip.

#### [NOTE]

There is no need to treat the adhering surface if it is non-precious metal.

#### B-3. Application of ED PRIMER II to the root canal and the cavity wall

For the application of ED PRIMER II to the root canal and the cavity wall, refer to A-5 in this chapter.

#### B-4. Preparing PANAVIA F 2.0 Paste

Refer to the foregoing section "Preparing PANAVIA F 2.0 Paste " in this chapter.

#### B-5. Placing the post

- Apply the mixed PANAVIA F 2.0 Paste to the post.
- After applying the mixed paste to the post, insert it into the root canal quickly. It is advisable to lightly vibrate the post while inserting it into the root canal to prevent the entrapment of air bubbles. Cementation should be completed within 40 seconds.

#### [CAUTION]

- Do not use a lentulo spiral for loading the PANAVIA F 2.0 Paste into the root canal. If the paste is loaded into the root canal using a composite resin syringe, the polymerization of the paste is accelerated. It is necessary to fit the post as quickly as possible.
- If multiple posts are to be fitted to a single tooth, use caution to prevent excess PANAVIA F 2.0 Paste from flowing into other root canals.

#### B-6. Spreading the excess PANAVIA F 2.0 Paste

Using a disposable brush tip, spread the excess paste over the coronal base and post head.

#### B-7. Light-curing

Light-cure PANAVIA F 2.0 Paste on the remaining crown and post head. See the tables "Dental curing unit" and "Light-curing time in A-9.

#### [CAUTION]

For the application of Opaque paste, light-cure is not necessary. Please proceed to the next step.

#### B-8. Core build-up

After cementing the post in the canal, place a core build-up composite resin (e.g. CLEARFIL PHOTO CORE or CLEARFIL DC CORE AUTOMIX) for preparation of the abutment tooth according to the manufacturer's instructions.

#### C. Standard procedure III (indications [6])

- Amalgam bonding

#### C-1. Cleaning of tooth structure

Clean the cavity and provide moisture control in the usual manner.

#### C-2. Application of ED PRIMER II to the cavity

For the application of ED PRIMER II to the cavity, refer to A-5 in this chapter.

#### C-3. Preparing PANAVIA F 2.0 Paste

Refer to the foregoing section "Preparing PANAVIA F 2.0 Paste " in this chapter.

#### C-4. Placing the amalgam

- Apply the mixed PANAVIA F 2.0 Paste to the cavity. Apply a thin, even layer of the mixed paste to the entire cavity surface primed with ED PRIMER II being careful to avoid entrapping air.

#### [CAUTION]

Because ED PRIMER II accelerates the set of PANAVIA F 2.0 Paste, the paste should be applied to the primed cavity quickly.

#### 2. Amalgam filling

The triturated amalgam should be condensed on the unset PANAVIA F 2.0 Paste. Occlusal carving can be accomplished in the normal manner.

#### C-5. Removing the excess PANAVIA F 2.0 Paste

Any excess paste remaining at the margins can be removed with a disposable brush tip.

#### C-6. Curing PANAVIA F 2.0 Paste

For curing the mixed paste along the cement margin, refer to A-9 in this chapter.

#### C-7. Finishing

Remove the excess PANAVIA F 2.0 Paste adhered to tooth structure by polishing.

#### [CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

#### [WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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"CLEARFIL", "PANAVIA", "OXYGUARD", "TRI-S BOND", "CLEARFIL PHOTO CORE" and "CLEARFIL DC CORE AUTOMIX" are trademarks of KURARAY CO., LTD.