kuraray



# CLEARFIL'SE Protect

**SE Protects your patients** 



### **Advanced New Features for More Confidence**


## **Antibacterial Cavity Cleansing Effect from MDPB**

Studies show that the risk of bacteria remaining in a cavity tends to increase with the smaller, "minimal intervention" cavities<sup>1)</sup>. Many doctors are using cavity cleansers, such as Consepsis (Ultradent) to address this tendency. CLEARFIL™ SE Protect contains a new functional monomer MDPB, which exhibits an "Antibacterial Cavity Cleansing Effect."

The antimicrobial effect with CLEARFIL™ SE Protect is greater than with current cavity cleansers as shown in Fig. 1, and is biologically safe because MDPB is fixed and will not leak after curing². In addition, CLEARFIL™ SE Protect is a cost and time saver by eliminating extra steps.

1) S. Imazato; Dent. Mater. J., 28(1), 11-19, 2009

2) E.A. Koulaouzidou et al.: J. Contep Dent. Pract., 9(7), 57-64, 2008



**Fig. 1** Number of Bacterial Recovery: Data from M. Türkün et al. Am. J. Dent., Vol. 19, 166-170, 2006.

\*Not a registered trademark of Kuraray Co., Ltd.




Fig. 2 The bactericidal mechanism of MDPB is presumed to be similar to the well-know antibacterial agent CPC (Cetyl pyridinium chloride), which is in many toothpastes and mouth rinses.

# Composite Adhesive **Test Specimen** Exposed interface was Demineralized with pH 4.5 solution for 90 Min SEM Observation Demineralized with pH 4.5 Acid-Base I solution for Resistant 90min. Zone (Super Dentin) A thickening due **Bond** to fluoride **Dentin** U X3,500

# Fluoride-Releasing

CLEARFIL™ SE Protect contains a patented™, specially treated sodium fluoride (NaF); the NaF in CLEARFIL™ SE Protect is coated with a unique polymer capsule that allows release of the NaF while the bonding layer physical properties, including strength, are maintained.

CLEARFIL™ SE Protect makes an "acid-base resistant zone" (called Super Dentin) interface between the bond and dentin behind the hybrid layer (Fig. 3, SEM image), as a result of the fluoride-release. Etch & rinse (total-etch) products like Adper Single Bond\* cannot make this zone<sup>3) 4)</sup>. This zone is 1.2µm which is thicker than the hybrid layer and very important for resisting the acid generated by bacteria if microleakage does occur<sup>3) 4)</sup>.

**Fig. 3** SEM image after soaking in a mild acid condition: Courtesy of G. Inoue, phases. T. Nikaido and J. Tagami, Tokyo Medical and Dental University, Japan

3) K. Waidyasekera, J. Tagami et al., J. Dent (37), 604-609, 2009 4) T. Nikaido, J. Tagami et al., Biomed Mater Eng. 2009;19(2-3):163-71

## Features from CLEARFIL™ SE BOND for Predictable Treatment

## **Virtually No Post-Operative Sensitivity**

Studies evaluating post-operative sensitivity of CLEARFIL™ SE Protect were published by Dr. Türkün, Ege Univ. School of Dentistry (Turkey) and THE DENTAL ADVISOR. The results reported that no postoperative sensitivity was experienced during one year<sup>5)</sup> and during 6 weeks<sup>6)</sup> evaluation period, respectively. The image in Fig. 4 was created using Transmission Electron Microscopy (TEM) technology. It clearly shows that CLEARFIL™ SE Protect bonds strongly to dentin. No gaps is a key factor in preventing post-operative sensitivity.

5) L. S. Türkün et al. JADA Vol. 136, 657-664 2005 6) The Dental Advisor, Vol. 21, No. 8, Oct. 2004





Fig. 4 TEM Micrograph: Courtesy of Dr. Hidehiko Sano, et al., Hokkaido Univ. Japan

# **Superior Bond Strength and Durability**

The test results of tensile bond strength (Fig. 5) show that CLEARFIL™ SE Protect has clear advantages regarding bond strength and durability compared with 3- or 2-step total-etch (etch and rinse) bonding agents, such as Adper Scotchbond Multi-Purpose\* (3-step) or Adper Single Bond\* (2-step), and 1-step self-etch such as i Bond\*. In addition, CLEARFIL™ SE Protect showed excellent clinical performance in high stress bearing areas for at least five years<sup>7)</sup>.



# **SE Protects Your Patients!**

#### Features from CLEARFIL™ SE BOND

- Virtually no post-operative sensitivity
- Excellent bond strength for enamel and dentin
- Fast and simple procedure with low technique sensitivity

#### **Advanced New Features**


- Antibacterial cavity cleansing effect
- Fluoride-releasing properties
- Higher, more durable bond than CLEARFIL™ SE BOND



# **Basic Procedure**


#### **Prepared Tooth**





#### **Primer**






- Apply PRIMER and leave for 20 sec.
- Dry with gentle air flow

#### **Bond**





- Apply BOND
- Air flow gently
- Light-cure for 10 sec.

#### Restoration



# Ordering information for CLEARFIL™ SE Protect



#### #2870-KA Standard Kit

1-Bond (6ml) 1-Primer (5ml)

#### Kit Accessories:

2-Brush tip handles 1-Mixing Dish 1-Light-Blocking Plate 1 Outer case

100-Disposable brush tips



#### #2872-KA Value Kit

3-Bond (6ml) 3-Primer (5ml)

#### Kit Accessories:

2-Brush tip handles 1-Mixing Dish

1-Light-Blocking Plate

1 Outer case

300-Disposable brush tips



#2881-KA Bond Bottle 1-Bond (5ml each)



#2882-KA
Primer Bottle
1-Primer
(6ml each)

#### YOUR CONTACT

#### Kuraray America, Inc.

32 Old Slip Phone (800) 879-1676
Floor 7 Fax (888) 700-5200
New York, NY 10005 E-mail info@kuraraydental.com
Website www.kuraraydental.com

- Before using this product, be sure to read the Instructions for Use supplied with the product.
- The specifications and appearance of the product are subject to change without notice.
- Printed color can be slightly different from actual color.

"CLEARFIL" is registered trademark or trademarks of KURARAY CO., LTD.





