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(TI-010)

## Dental Zirconia

## HT (High Translucent) & ML (Multi Layered)

### ENGLISH INSTRUCTIONS FOR USE

#### I. Introduction


KATANA Zirconia has 2 shapes available: Block and Disc. This IFU is for KATANA Zirconia HT (High Translucent) & ML (Multi Layered). KATANA Zirconia HT & ML are a pre-sintered zirconia disc 98.5 mm in diameter, which contains a plastic ring. This is designed for all milling systems using this generic-type disc. (Please refer to your milling system’s technical instructions for correct machine operation.) KATANA Zirconia HT has five available thicknesses: (10mm, 14mm, 18mm, 22mm and 26mm) and has three shade variations: (HT10, HT12 and HT13). KATANA Zirconia ML has three available thicknesses: (14mm, 18mm and 22mm) and has six shade variations: (A Light, A Dark, B Light, C Light, D Light and A White). KATANA Zirconia HT is higher in translucency and KATANA Zirconia ML consists of four gradated shade layers. KATANA Zirconia HT & ML are recommended for use in fabricating FCZ (Full Contour Zirconia) restorations or the frameworks.

#### II. Intended Use

KATANA Zirconia is used for the fabrication of the all-ceramic restorations (frameworks, FCZ crowns, FCZ bridges, inlays, onlays and veneers).

#### III. Sintering Program

Temperature	Programming Rate	Holding Time
Room Temp. — 1500°C (2732°F)	10°C/min. (18°F/min.)	—
1500°C (2732°F)	—	2 hrs
1500°C (2732°F) — Room Temp.	-10°C/min. (-18°F/min.)	—

#### IV. Composition

ZrO<sub>2</sub>, Y<sub>2</sub>O<sub>3</sub> etc.

#### V. Type and Class (ISO6872:2008)

Type:II/ Classe:6

#### VI. Physical Properties

Coefficient of Thermal Expansion (25-500°C (77-932°F)): 9.9(±0.2)×10<sup>-6</sup>/K

#### VII. Directions for Use

- Take the disc from the packaging and confirm that the disc does not have a crack or other damage.
- Place the disc into the milling machine; then begin the milling process following the milling systems technical instructions.
- After milling, remove the restorations from the disc with a diamond bur etc.

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
- The cutting waste or dust, which is attached to the restorations, can be removed with a gentle air stream.
- Put the restorations into the refractory sagger tray and place them in the sintering furnace.
- Sinter the restorations according to the sintering schedule shown above. (III. Sintering Program).
- After sintering, adjust the restorations with a diamond bur as needed.
- Confirm that the restorations have no cracks.
- 1 FCZ:
  - Glaze baking: Create a high shine surface by polishing, especially on the contact areas then apply the glaze on all surfaces in the usual manner.
  - Hand polishing: Create a high shine surface on the entire restoration by polishing without using the glaze.
- 2 Frameworks: Build-up the dental porcelain (CERABIEN ZR o CZR PRESS LF) or press the ingot (CZR PRESS) on the frameworks following the manufacturer’s technical instructions.

#### VIII. Remarks on Handling

##### Contraindicaciones:

• If the patient is hypersensitive to zirconia or any other components, this product must not be used.

##### Warning:

• If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis etc., discontinue use of the product and seek medical attention immediately.

##### Caution:

- When milling the disc, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.
- When milling the disc, use safety glasses to prevent the dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.
- Do not use for any purposes except for dental restoration. This product is for dental application only.
- Do not touch the items heated by the furnace with your bare hands.

##### Caution in conjunction with usage:

- Do not use the disc if there is a crack noticed after removing it from the package.
- Dry milling is recommended. If wet milling is used, the translucency of the disc may be reduced.
- If you find a crack in the restorations, do not use.
- The above sintering recommendation is only a guideline; some adjustments may be required depending on each individual furnace.
- When milling the disc, do not bring the milling bar into contact with the plastic ring to prevent detachment of the disc.
- Cut and remove the plastic ring prior to sintering the whole disc.

##### Storage:

- Store in a cool and dry place. Keep away from direct sunlight.
- The product should be stored at 10-30°C (50-86°F).
- Do not remove the disc from its packaging during storage.
- The disc is fragile, and requires care when handling.
- The product must be stored in an appropriate place where only dental personnel have access.
- The product must be used by the expiration date indicated on the package.

#### [WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

#### [NOTE]

"KATANA" and "CERABIEN" are trademarks of NORITAKE CO., LIMITED.

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Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a dental professional.
Kuraray Europe GmbH