# CLEARFIL™ CERAMIC PRIMER

## ENGLISH INSTRUCTIONS FOR USE

### I. INTRODUCTION

CLEARFIL CERAMIC PRIMER is a silane-coupling agent that provides an enhanced adhesive surface to porcelain, ceramics, hybrid ceramics or composite resin.

### II. INDICATIONS

CLEARFIL CERAMIC PRIMER is indicated for the following uses:

- [1] Surface treatment of porcelain, ceramics, hybrid ceramics or composite resin
- [2] Intraoral repairs of fractured crowns/bridges made of porcelain, ceramics, hybrid ceramics or composite resin

### III. CONTRAINDICATION

Patients with a history of hypersensitivity to methacrylate monomers.

### IV. POSSIBLE SIDE EFFECT

The oral mucosa membrane may turn whitish when contacted by the product due to the coagulation of protein. This is a temporary phenomenon that will usually disappear in a few days.

### V. PRECAUTIONS

#### 1. Safety precautions

1. Avoid using the product for patients with a history of hypersensitivity to methacrylate monomers.
2. If any hypersensitivity occurs, such as dermatitis, discontinue the use of the product and consult a physician.
3. Use caution to prevent the product from coming in contact with the skin or getting in the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from any splashing material.
4. If the product comes in contact with the human body, take the following actions:
  - <If the product gets in the eye>  
Immediately wash the eye with copious amounts of water and consult a physician.
  - <If the product comes in contact with the skin and the oral mucosa>  
Immediately wipe it off with a cotton pledget or gauze moistened with alcohol and then wash with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. Avoid using the same disposable brush tip for different patients to prevent cross contamination. Discard the tip after use and sterilize the brush tip handle.
7. Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers.

#### 2. Handling and manipulation precautions

1. Do not mix the product with bonding agents or use it in conjunction with other silane-coupling agent; intermixing materials may cause a change in physical properties, including a possible decrease from the properties expected.
2. Do not use the product near an open flame. It contains ethanol, a flammable substance.
3. The product should be used as soon as practicable after dispensing. The product contains volatile ethanol. As the solvent evaporates, the viscosity increases and could make it difficult to apply.
4. Perform bonding (or cementation) soon after treating the restoration surface with the product.
5. If the treated surface is contaminated, wash it with water, dry, clean with an etching agent (e.g. K-ETCHANT GEL) treatment according to the manufacturer's instructions, and treat again with the product.
6. If the product gets deposited on the nozzle or the container, wipe it away using a piece of gauze or a cotton pad soaked with alcohol.
7. Use caution not to get your fingers injured by any sharp edge that might occur due to damage.
8. The use of the product is restricted to a licensed dentist.
9. The product should not be used for purposes other than specified in [II. INDICATIONS].

#### 3. Storage precautions

1. Do not use the product after the expiration date. Note the expiration date on the outside of the package.
2. The product should be refrigerated (2-8°C/36-46°F) when not in use, and should be brought to room temperature before using. After taking out of the refrigerator, the product must be left standing until it comes to room temperature; otherwise, the liquid might drop more than necessary or ooze after use.
3. Keep away from extreme heat, direct sunlight, or flames.
4. Replace the bottle cap as soon as possible after the product has been dispensed from the bottle. The volatile component will evaporate if the bottle is left uncapped.
5. The product must be stored in a proper location to assure that no one other than dental practitioners can gain access and use the product.

### VI. COMPONENTS

#### CLEARFIL CERAMIC PRIMER

Principal ingredients:

- 3-Methacryloxypropyl trimethoxy silane
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- Ethanol

### VII. CLINICAL PROCEDURES

#### A. Surface treatment of porcelain, ceramics, hybrid ceramics or composite resin

##### A-1. Moisture control

If necessary, place a rubber dam to accomplish an appropriate isolation for adhesive procedures.

##### A-2. Surface preparation

If necessary, sandblast the adherent surface using alumina particles. After sandblasting, clean the adherent surface in an ultrasonic unit for 2 minutes.

##### A-3. Phosphoric acid treatment

Apply an etching agent (e.g. K-ETCHANT GEL) to the adherent surface according to the manufacturer's instructions and leave it in place for 5 seconds. The surface is then washed with water and dried. There is no need to treat if the adherent surface is ceramic consisting primarily of zirconia or alumina.

##### A-4. Application of CLEARFIL CERAMIC PRIMER

Dispense the necessary amount of CLEARFIL CERAMIC PRIMER into a well of the mixing dish immediately before application. Apply CLEARFIL CERAMIC PRIMER to the adherent surface with a disposable brush tip. After application, dry the entire adherent surface sufficiently by blowing mild oil-free air.

##### A-5. Cementation

Cement the prosthetic appliance using resin cement (e.g. PANAVIA F 2.0) according to the manufacturer's instructions.

#### B. Intraoral repairs of fractured crowns/bridges made of porcelain, ceramics, hybrid ceramics or composite resin

##### [NOTE]

The following are the clinical procedures for using the Clearfil tri-S bond Single dose system. If other bond systems are to be used, follow the Instructions for Use supplied with the bond system you are going to use.

##### B-1. Moisture control

If necessary, place a rubber dam to accomplish an appropriate isolation for adhesive procedures.

##### B-2. Preparation of fractured surfaces

1. Facing material surface  
Using a diamond point, remove a layer of the fractured surface and roughen the adherent surface. If necessary, place a bevel at the marginal area.
2. Exposed surface of metal frame  
Roughen the metal surface with a diamond point.

##### B-3. Acid etching of facing material and exposed metal surface

Apply an etching agent (e.g. K-ETCHANT GEL) to the adherent surfaces of the facing material and the metal frame according to the manufacturer's instructions and leave it in place for 5 seconds. The surfaces are then washed with water and dried.

##### B-4. Treatment of exposed metal surface

Apply a metal-adhesive primer (e.g. ALLOY PRIMER) to the metal surface according to the manufacturer's instructions.

##### B-5. Silane treatment of facing material surface

Dispense the necessary amount of CLEARFIL CERAMIC PRIMER into a well of the mixing dish immediately before application. Apply CLEARFIL CERAMIC PRIMER to the facing material surface with a disposable brush tip. After application, dry the entire adherent surface sufficiently by blowing mild oil-free air.

##### B-6. Treatment of tooth surface

When there is tooth surface as part of the adherent surfaces, apply a bonding agent (e.g. Clearfil tri-S bond Single dose) to it according to the manufacturer's instructions.

##### B-7. Light-cure composite resin filling

Apply a light-cure composite resin (e.g. CLEARFIL AP-X) into the cavity, light-cure, finish and polish the resin according to the manufacturer's instructions.

##### [CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

##### [WARRANTY]

KURARAY MEDICAL INC. will replace any product that is proved to be defective. KURARAY MEDICAL INC. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

##### [NOTE]

CLEARFIL and PANAVIA are trademarks of KURARAY CO., LTD.

Manufactured by  
**KURARAY MEDICAL INC.**  
1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

Distributed by  
**KURARAY AMERICA, INC.**  
600 Lexington Avenue, 26th Floor, New York NY 10022  
(800)879-1676

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CLEARFIL y PANAVIA son marcas registradas de KURARAY CO., LTD.