# CLEARFIL™ DC BOND

## ENGLISH INSTRUCTIONS FOR USE

### I. INTRODUCTION

CLEARFIL DC BOND is a dual-cure (light- and/or self-cure), two-component, one-step bonding system that allows simultaneous treatment of both dentin and enamel.

### II. INDICATIONS

- CLEARFIL DC BOND is indicated for the following uses:
- Core build-ups using light-, dual-, or self-cure composite resin
  - Direct restorations using light-cure composite resin
  - Direct restorations using self-cure composite resin
  - Cavity sealing as a pretreatment for indirect restorations
  - Treatment of exposed root surfaces
  - Intraoral repairs of fractured crowns/bridges made of ceramics, hybrid ceramics or composite resin

### III. CONTRAINDICATION

Patients with a history of hypersensitivity to methacrylate monomers

### IV. POSSIBLE SIDE EFFECT

The oral mucosa membrane may turn whitish when contacted by the product due to the coagulation of protein. This is a temporary occurrence that will usually disappear in a few days.

### V. INCOMPATIBILITIES

- Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol could retard the bonding system curing process.
- Do not use hemostatic agents, especially those containing ferric compounds, since these materials may impair adhesion and may cause discoloration at the tooth margin or surrounding gingiva due to ferric ions that may remain.

### VI. PRECAUTIONS

- Safety precautions**
  - Avoid using the product on patients with a known history of hypersensitivity to methacrylate monomers.
  - If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching, or numbness, discontinue the use of the product and seek medical attention.
  - Use caution when using the product to prevent contact with the soft oral tissue or skin. If the product comes in contact with the soft oral tissue or skin, wipe away with a cotton pledget moistened with alcohol and immediately rinse with copious amounts of water. Prior to use, cover the patient's eyes with a towel or safety glasses to protect from splashing materials. If the product gets in the eye, immediately rinse with copious amounts of water and consult an ophthalmologist.
  - Use caution when using the product to prevent swallowing.
  - Avoid looking directly at the curing light; take necessary protective measures.
  - Do not use the same Disposable brush tip for different patients to prevent cross contamination. Discard the tip after use and sterilize the brush tip handle after each patient.
  - Avoid direct contact with the skin and / or soft tissue to prevent hypersensitivity. Wear gloves or take appropriate precautions when using the product.

### 2. Handling and manipulation precautions

- Do not use a spiral filler (tenuilo spiral) for loading the paste into the root canal.
- Dispense BOND Liquid A and B from the containers by pointing the container downward and as vertically as possible; careful dispensing is necessary otherwise the liquid amounts may be unequal and the function of BOND will be impaired.
- When dispensing BOND Liquid B, press the container gently; otherwise the liquid might be dispensed in large quantities since the viscosity of BOND Liquid B is lower than BOND Liquid A.
- BOND Liquid A and B should be mixed when in use. Do not use them separately as a single agent.
- Do not use the product in conjunction with other bonding agents. Intermingling materials may cause a change in physical properties, including a possible decrease from the properties expected.
- Use Light blocking plate to avoid exposing BOND to an operating light or natural light. In addition, during tooth treatment with BOND for 20 seconds, move the lighting-spot out of the mouth or turn off the light to prevent the applied BOND from being exposed to strong light. BOND will set to a gel if left under an operating light or natural light (sunlight from windows).
- Use the mixture of BOND Liquid A and B as soon as possible after mixing. BOND contains volatile ethanol that as the solvent evaporates, the viscosity increases making it difficult to apply. The mixture must be covered with Light blocking plate and used within 90 seconds after mixing.
- Do not use BOND near an open flame. It contains ethanol, a flammable substance.
- Any overfilled resin that remains on uncut enamel should be removed thoroughly; marginal discoloration could result from the overfilled resin tags. In such case that there is a possibility of resin spreading over uncut enamel but the removal seems difficult, apply an etching agent (e.g. K-ETCHANT GEL) to the enamel according to the manufacturer's instructions, let it remain for 10 seconds, wash with water, and then dry.
- Dry the entire adherent surface sufficiently using air from a air syringe; the adhesion effect could be impaired if the surface is not dried sufficiently.
- Do not use the product for the surface treatment of metal posts.
- The emitting tip of the dental curing light should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light-evaluating device at appropriate intervals. Low light intensity causes poor adhesion.
- Be careful not to cut your fingers on the sharp edges of the instruments.
- The use of the product is restricted to a licensed dentist.
- Do not use the product for any purposes other than those specified in INDICATIONS.

### 3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- The product must be refrigerated (2-8°C / 36-46°F) when not in use, and should be brought to room temperature before using. After taking out of the refrigerator, the product must be left standing for more than 15 minutes, or until it comes to room temperature; otherwise bubbles will form in the liquid when dispensing or ooze after use.
- Keep away from extreme heat, direct sunlight or flame.
- Cap the BOND container immediately and completely after dispensing the liquid. If the container is left with the cap improperly replaced, the volatile substances will evaporate and the function of BOND will be impaired. If the liquid does not flow easily from the nozzle it may have

- plugged up the nozzle valve. Do not try to dispense the liquid forcibly.
- Wipe any BOND Liquid B off the nozzle of its container immediately after dispensing the liquid; otherwise white or pale yellow crystals will get deposited on the nozzle. If crystals get deposited on the nozzle, wipe them off cotton or gauze moistened with ethanol.
- The product must be carefully stored and used by a licensed dental professional.

### VII. COMPONENTS

Please see the outside of the package for quantity.

#### 1) BOND Liquid A

- Principal ingredients:
- 2-Hydroxyethyl methacrylate (HEMA)
  - Bis-phenol A diglycidylmethacrylate (Bis-GMA)
  - 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
  - di-Camphorquinone
  - Benzoyl peroxide
  - Colloidal silica

#### 2) BOND Liquid B

- Principal ingredients:
- Water
  - Ethanol

#### 3) Accessories

- Disposable brush tips
- Brush tip handle
- Mixing dish
- Light blocking plate

### VIII. CLINICAL PROCEDURES

#### A. Core build-ups using light-, dual-, or self-cure composite resin

##### 1-1. Cleaning tooth structure

An adequately cleaned tooth surface assures maximum adhesive performance. Be sure the tooth surface is adequately cleaned.

##### A-2. Moisture control

In order to produce optimal results, avoid contamination of the treatment area from saliva or exudates. A rubber dam is recommended to keep the tooth clean and dry.

##### A-3. Tooth and root canal preparation

Remove existing restorations, decay and/or caries from the tooth and prepare root canal in the usual manner.

##### A-4. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

##### A-5. Post preparation

Select an appropriate post and make necessary adjustments and preparation to it.

##### A-6. Treatment of post surfaces (Either A-6a or A-6b)

- A-6a. For metal posts**
  - If the dental post is not already surface treated, sandblast it with alumina powder.
  - Apply a metal-adhesive primer (e.g. ALLOY PRIMER) to the post surface with a brush, according to the manufacturer's instructions.

##### A-6b. For glass-fiber posts

- Apply an etching agent (e.g. K-ETCHANT GEL) to the post surface according to the manufacturer's instructions and leave it in place for 5 seconds. Then wash the surface with water and dry.
- Apply a silane-coupling agent (e.g. CLEARFIL CERAMIC PRIMER) to the post surface according to the manufacturer's instructions.

#### [CAUTION]

Do not use sand blast alumina with glass-fiber posts; this could damage the post fibers. Damage to the fibers degrades the bond strength to adhesive resin materials.

#### A-7. Application of BOND

- Dispense equal amounts of BOND Liquid A and B into a well of Mixing dish and mix for more than 5 seconds immediately before application.

#### [CAUTION]

Use the mixture of BOND Liquid A and B as soon as possible after mixing. The mixture must be covered with Light blocking plate and used within 90 seconds after mixing.

- Apply mixed BOND to the root canal and the cavity wall with a Disposable brush tip. Leave it in place for 20 seconds. Use caution not to allow saliva or exudates to contact the treated surfaces for at least 20 seconds.

#### [CAUTION]

BOND will set to a gel if left under an operating light. Move the lighting-spot out of the mouth or turn off the light to prevent the applied BOND from being exposed to strong light.

- After conditioning the adherent surface for 20 seconds, dry the entire adherent surface sufficiently by blowing high-pressure, oil-free air for more than 5 seconds while spreading the bond layer thinly. Use a vacuum aspirator to prevent the BOND liquid from scattering. Remove excess BOND with a paper point and once again dry the adherent surface sufficiently by blowing high-pressure, oil-free air.

#### [CAUTION]

BOND contains water and volatile ethanol. Dry the entire adherent surface sufficiently by blowing oil-free air; otherwise optimal adhesion will be impaired. Observe the drying method and treatment time to ensure optimum adhesion.

- If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat with BOND again.

- Light-cure BOND with a dental curing light (see table "Dental curing light") for the specified length of time shown in table "Light-curing time for dental curing lights".

Table: Dental curing light

Type	Light source	Wavelength range and light intensity
Conventional halogen	Halogen lamp	Light intensity <sup>2)</sup> of 300 - 550 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm
Fast halogen	Halogen lamp	Light intensity <sup>2)</sup> of more than 550 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm
Plasma arc	Xenon lamp	Light intensity <sup>3)</sup> of more than 2000 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm, and light intensity of more than 450 mW/cm <sup>2</sup> in wavelength range from 400 - 430 nm.
LED	Blue LED <sup>1)</sup>	Light intensity <sup>3)</sup> of more than 300 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm

1) Peak of emission spectrum: 450 - 480 nm

2) Evaluated according to ISO 10651-1

3) Wavelength distribution and light intensity values measured with a spectro-radiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp

Table: Light-curing time for dental curing lights

Dental curing light	Light-curing time
Conventional halogen	20 seconds
LED	20 seconds
Fast halogen	10 seconds
Plasma arc	10 seconds

#### A-8. Post cementation

- Apply composite resin for core build-up (e.g. CLEARFIL DC CORE AUTOMIX) into the root canal according to the manufacturer's instructions.
- Insert the post into the root canal and fix it in place.

#### [CAUTION]

- When using a dual-cure composite resin for core build-up, light-cure using the correct dental curing light, for the specified length of time, when fixing the post in place. Be sure to light-cure the margins of the tooth structure and the filling composite resin sufficiently to assure a good bond strength.
- When using a self-cure composite resin for core build-up, be sure to leave it for the specified length of time for complete hardening.
- The chemical polymerization of composite resin for core build-up (e.g. CLEARFIL DC CORE AUTOMIX) is accelerated on contact with a surface on which the BOND has been applied.

#### A-9. Core build-up

- After cementing the post in the root canal, place the core build-up composite resin around the post.
- Complete the curing process according to the manufacture's instructions.

#### [CAUTION]

When using a light-cure composite resin for core build-up, be sure to light-cure the material. Check the polymerization depth of the composite resin in the Instructions for Use and apply the composite resin in a thickness that allows the light to reach the bond layer through the paste. After building up the composite resin, light-cure it from both the lingual and the labial sides to assure optimal bond strength.

When using a dual-cure composite resin for core build-up, check the polymerization depth of the composite resin in the Instructions for Use and apply the composite resin in a thickness that allows the light to reach the bond layer through paste. Light-cure it from both the lingual and the labial sides to assure optimal bond strength. Considering the thickness, leave it for the specified length of time for additional self-curing after the light curing.

- When using a self-cure composite resin for core build-up, place the composite resin and leave it for the specified curing time.
- A-10. Preparing an abutment tooth**  
After being sure that the core build-up composite resin is completely cured, finish the construction of the core in the usual manner.

**B. Direct restorations using light- or self-cure composite resin / cavity sealing and treatment of exposed root surfaces**

##### B-1. Cleaning tooth structure

An adequately cleaned cavity assures maximum adhesive performance. Be sure the cavity is adequately cleaned.

##### B-2. Moisture control

In order to produce optimal results, avoid contamination of the treatment area from saliva or exudates. A rubber dam is recommended to keep the tooth clean and dry.

##### B-3. Cavity preparations

Remove any infected dentin using a caries detection material (e.g. CARIES DETECTOR) according to the manufacturer's instructions, and prepare the cavity in the usual manner.

##### B-4. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

##### B-5. Acid etching uncut enamel

If there is a possibility of resin spreading over uncut enamel, apply an etching agent (e.g. K-ETCHANT GEL) to the uncut enamel according to the manufacturer's instructions, let it remain for 10 seconds, wash with water, and then dry.

##### B-6. Application of BOND

- Dispense equal amounts of BOND Liquid A and B into a well of Mixing dish and mix for more than 5 seconds immediately before application.

#### [CAUTION]

Use the mixture of BOND Liquid A and B as soon as possible after mixing. The mixture must be covered with Light blocking plate and used within 90 seconds after mixing.

- Apply mixed BOND to the entire cavity wall with a Disposable brush tip. Leave it in place for 20 seconds. Use caution not to allow saliva or exudates to contact the treated surfaces for at least 20 seconds.

#### [CAUTION]

BOND will set to a gel if left under an operating light. Move the lighting-spot out of the mouth or turn off the light to prevent the applied BOND from being exposed to strong light.

- After conditioning the adherent surface for 20 seconds, dry the entire adherent surface sufficiently by blowing high-pressure oil-free air for more than 5 seconds while spreading the bond layer thinly. Use a vacuum aspirator to prevent the bond liquid from scattering.

#### [CAUTION]

BOND contains water and volatile ethanol. Dry the entire adherent surface sufficiently by blowing oil-free air; otherwise optimal adhesion will be impaired. Observe the drying method and treatment time to ensure optimum adhesion.

- If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat with BOND again.

- Light-cure BOND with a dental curing light (see table "Dental curing light") for the specified length of time shown below.

Table: Light-curing time for dental curing lights

Dental curing light	Light-curing time
Conventional halogen	10 seconds
LED	10 seconds
Fast halogen	5 seconds
Plasma arc	5 seconds

#### B-7. Follow either B-7a, B-7b or B-7c

##### B-7a. Direct restorations using light-cure composite resin

Apply composite resin (e.g. CLEARFIL MAJESTY Esthetic, CLEARFIL MAJESTY Posterior) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

##### B-7b. Direct restorations using self-cure composite resin

Apply composite resin (e.g. CLEARFIL F II NEW BOND) into the cavity, finish and polish according to the manufacturer's instructions.

##### B-7c. Cavity sealing and treatment of exposed root surfaces

Apply a thin coat of flowable light-cure composite resin (e.g. CLEARFIL MAJESTY Flow) to the tooth, and light-cure it according to the manufacturer's instructions. Remove unpolymerized resin with a cotton pledget moistened with alcohol.

### C. Intraoral repairs of fractured crowns/bridges made of ceramics, hybrid ceramics composite resin

#### C-1. Preparation of fractured surface

- Facing material surface
  - Using a diamond point, remove a layer of the fractured surface and roughen the adherent surfaces. If necessary, place a bevel at the marginal area.
  - Exposed surface of metal frame
    - Roughen the metal surface with a diamond point.

#### C-2. Acid etching of facing material and metal surfaces

Apply an etching agent (e.g. K-ETCHANT GEL) to the adherent surfaces of the facing material and the metal frame according to the manufacturer's instructions and leave it in place for 5 seconds. Then wash the surfaces with water and dry.

#### C-3. Treatment of exposed metal surface

Apply a metal-adhesive primer (e.g. ALLOY PRIMER) to the metal surface according to the manufacturer's instructions.

#### C-4. Silane treatment of facing material surface

Apply a silane-coupling agent (e.g. CLEARFIL CERAMIC PRIMER) to the facing material surface with Disposable brush tip, according to the manufacturer's instructions.

#### C-5. Application of BOND (When there is tooth surface as par of the adherent surface)

- Dispense equal amounts of BOND Liquid A and B into a well of Mixing dish and mix for more than 5 seconds immediately before application.

#### [CAUTION]

Use the mixture of BOND Liquid A and B as soon as possible after mixing. The mixture must be covered with Light blocking plate and used within 90 seconds after mixing.

- Apply the mixed BOND to the adherent surface with a Disposable brush tip. Leave it in place for 20 seconds. Use caution not to allow saliva or exudates to contact the treated surfaces for at least 20 seconds.

#### [CAUTION]

BOND will set to a gel if left under an operating light. Move the lighting-spot out of the mouth or turn off the light to prevent the applied BOND from being exposed to strong light.

- After conditioning the adherent surface for 20 seconds, dry the entire adherent surface sufficiently by blowing high-pressure oil-free air for more than 5 seconds while spreading the bond layer thinly. Use a vacuum aspirator to prevent the bond liquid from scattering.

#### [CAUTION]

BOND contains water and volatile ethanol. Dry the entire adherent surface sufficiently by blowing oil-free air; otherwise the adhesion effect will be impaired. Observe the drying method and treatment time to ensure optimal adhesion.

- If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat with BOND again.

- Light-cure BOND with a dental curing light (see table "Dental curing light" in A-7.), for the specified length of time shown below.

Table: Light-curing time for dental curing lights

Dental curing light	Light-curing time
Conventional halogen	10 seconds
LED	10 seconds
Fast halogen	5 seconds
Plasma arc	5 seconds

#### C-6. Light-cure composite resin filling

Apply light-cure composite resin (e.g. CLEARFIL MAJESTY Esthetic, CLEARFIL MAJESTY Posterior) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

#### [NOTE]

If necessary, use opaque resin such as CLEARFIL ST OPAQUER for the metal surface to prevent metal from shining through.

#### [CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

#### [WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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CLEARFIL, CLEARFIL DC CORE AUTOMIX, CLEARFIL MAJESTY and CLEARFIL ST are trademark of KURARAY CO., LTD.