CLEARFIL™ SE BOND Restorative Kit

ENGLISH INSTRUCTIONS FOR USE

I. INTRODUCTION

CLEARFIL AP-X PLT is a light-cured, radiopaque restorative composite resin which provides accurate colormatching, high polishability and excellent physical properties, making it ideal for both anterior and posterior restorations. It is formulated with optimal viscosity assuring easy handling and placement. CLEARFIL AP-X PLT, with the special dispensing system, can be quickly placed directly into the cavity. CLEARFIL SE BOND is a light-cured bonding system that consists of a self-etching primer and a bonding agent. The primer provides simultaneous treatment of both dentin and enamel.

II. INDICATIONS

CLEARFIL SE BOND Restorative Kit is indicated for the following restorative applications:

- Class I, II, V restorations of posterior teeth
- Class III, IV, V restorations of anterior teeth
- Cervical cavities or roots surfaces defects

III. CONTRAINDICATION

Patients with a history of hypersensitivity to methacrylate monomers

IV. POSSIBLE SIDE EFFECT

The oral mucosa membrane may turn whitish when contacted by the PRIMER due to the coagulation of protein. This is a temporary phenomenon that will usually disappear in a few days.

V. INCOMPATIBILITIES

- Do not use eugenol containing materials for pulp protection or temporary sealing, since the eugenol could retard the bonding system curing process.
- Do not use hemostatics containing ferric compounds, since these materials may impair adhesion and may cause discoloration in the tooth margin or surrounding gingival due to remaining ferric ions.

VI. PRECAUTIONS

1. Safety precautions

- Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomers.
- If any hypersensitivity occurs, such as dermatitis, discontinue use of the product and consult a physician.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
- If the product comes in contact with human body tissues, take the following actions:
 - < If the product gets in the eye > Immediately rinse the eye with copious amounts of water and consult a physician.
 - < If the product comes in contact with the skin or the oral mucosa > Immediately wipe it off with a cotton pledget moistened with alcohol or gauze and rinse with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Avoid looking directly at the dental curing light when curing the product.
- Avoid using the same Disposable brush tip for different patients to prevent cross infection. Discard the tip after use and sterilize the Brush tip handle.

2. Handling and manipulation precautions

- Do not use CLEARFIL AP-X PLT in conjunction with the other composite resin. Mixing materials may cause a change in physical properties, possibly a decrease, from the properties expected.
- When light curing CLEARFIL AP-X PLT, note the light curing depth in this Instructions For Use.
- CLEARFIL AP-X PLT and CLEARFIL SE BOND polymerize by visible light (especially ultraviolet rays). Avoid operating light or natural light (sunlight from windows) and use within 3 minutes of dispensing.
- The emitting tip of the dental curing light should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light-evaluating device at appropriate intervals.
- The use of the product is restricted to licensed dentist.

3. Storage precautions

- Do not use after the expiration date. Note the expiration date on the outside of package.
- CLEARFIL AP-X PLT must be stored below 25°C / 77°F when not in use. CLEARFIL SE BOND must be refrigerated (2-8°C / 36-46°F) when not in use, and should be brought to room temperature before using. After taking out of the refrigerator, the product must be left standing until it comes to room temperature; otherwise the liquid might drop more than necessary or ooze after use.
- Keep away from extreme heat, direct sunlight or flame.
- Replace the bottle cap as soon as possible after resin has been dispensed from the bottle.
- The product must be stored in proper places where only dental practitioners can access it.

VII. SYSTEM SHADES AND COMPONENTS

1. Shades

CLEARFIL AP-X PLT is available in 8 shades; 6 shades correspond on a best-match basis with VITA Shade Guide plus 2 shades for special requirements. Select an appropriate shade by matching the shade guide to the restoration site. The XL and HO shades offer a strong, lighter white color. The HO shade is transcendent white and is suitable for sites requiring extra white color. Standard shades: A2, A3, A3.5, B2, B3 Cervical shades: A4 Transcendent shades: XL, HO

2. Kit Components

CLEARFIL SE BOND Restorative Kit consists of CLEARFIL AP-X PLT, CLEARFIL SE BOND (PRIMER and BOND) and accessories.

1) CLEARFIL AP-X PLT
10tips (0.2g/tip) : A2, A3, B2, XL
Principal ingredients:

- Silanated barium glass powder
- Silanated colloidal silica
- Silanated α -quartz powder
- Bis-phenol A diglycidylmethacrylate (Bis-GMA)
- Triethylenglycol dimethacrylate (TEGDMA)
- di-Camphorquinone

[NOTE]

All 8 shades are available separately.

2) CLEARFIL SE BOND: PRIMER (6ml)

Principal ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- 2-Hydroxyethyl methacrylate (HEMA)
- di-Camphorquinone
- Water

3) CLEARFIL SE BOND: BOND (5ml)

Principal ingredients:

- Bis-phenol A diglycidylmethacrylate (Bis-GMA)
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- 2-Hydroxyethyl methacrylate (HEMA)
- di-Camphorquinone
- Silanated colloidal silica

4) Accessoires

- Outer case 1 piece
- PLT tip case 4 pieces

VIII. CLINICAL PROCEDURES

1. Cleaning tooth structure

Be sure the cavity is adequately cleaned. An adequately cleaned cavity assures maximum adhesive performance.

2. Moisture control

Avoid contamination of the treatment area from saliva or blood, in order to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

3. Cavity preparations

Remove any infected dentin using CARIES DETECTOR as a guide and prepare the cavity in the usual manner.

4. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

5. Shade selection

Before isolating the teeth, select an appropriate shade using the VITA Lumin VACUUM shade guide.

6. Acid etching uncut enamel

If there is a possibility of resin spreading over uncut enamel, apply K-ETCHANT GEL to the enamel, let it remain for 10 seconds, rinse with water, and then dry.

[CAUTION]

The use of PRIMER alone may not sufficiently

condition uncut enamel. Overfilling of resin onto unetched, uncut enamel could cause marginal discoloration.

7. Tooth surface treatment

- Dispense the necessary amount of PRIMER into a well of the mixing dish.
- Apply PRIMER to the entire cavity wall with a sponge or Disposable brush tip. Leave it in place for 20 seconds. Use caution not to allow saliva or exudates.
- After conditioning the tooth surface for 20 seconds, dry the entire adherent surface sufficiently by blowing mild oil-free air. Use a vacuum aspirator to prevent the PRIMER liquid from scattering.

[CAUTION]

- PRIMER will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing PRIMER to operating light or natural light, and use within 3 minutes from dispensing.
- Observe the drying method and treatment time to ensure optimum adhesion.
- Do not rinse treated surface with water, and avoid touching the treated surface. However if the treated surface is contaminated, rinse it with water, dry, or clean with alcohol, and treat with PRIMER again.

8. Bonding

- Dispense the necessary amount of BOND into a well of the mixing dish.
- Apply BOND to the entire cavity wall with a sponge or Disposable brush tip.
- After application, make the BOND film as uniform as possible using a gentle oil-free air stream.

[CAUTION]

- BOND will set to a gel if left under an operating light or natural light (sunlight from windows). Use the Light blocking plate to avoid exposing BOND to an operating light or natural light, and use within 3 minutes from dispensing.
- A strong air blast will scatter BOND, resulting in poor adhesion.

- Light-cure BOND with a dental curing unit* for 10 seconds (Conventional halogen, LED, Fast halogen, and Plasma arc).

* Dental curing unit

Type	Light source	Wavelength range and light intensity
Conventional halogen	Halogen lamp	Light intensity ²⁾ of 150 - 550 mW/cm ² in wavelength range from 400 - 515 nm
Fast halogen	Halogen lamp	Light intensity ²⁾ of more than 550 mW/cm ² in wavelength range from 400 - 515 nm
Plasma arc	Xenon lamp	Light intensity ²⁾ of more than 2000 mW/cm ² in wavelength range from 400 - 515 nm, and light intensity of more than 450 mW/cm ² in wavelength range from 400 - 430 nm.
LED	Blue LED ¹⁾	Light intensity ²⁾ of more than 120 mW/cm ² in wavelength range from 400 - 515 nm

¹⁾ Peak of emission spectrum: 450 - 480 nm

²⁾ Wavelength distribution and light intensity values measured with a spectro-radiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp

9. Placement and light curing of CLEARFIL AP-X PLT

- Dispensing
Place a PLT (Pre-loaded tip) into the dispenser barrel according to the instructions for use for dispenser. (Centrix syringe such as CR E/Z or MARK IIIp is recommended.)
The tip may be rotated to gain the proper angle for entrance into the cavity.
Squeeze the dispenser with a slow, steady pressure. Excessive force is not necessary. Remove and discard the tip after use.
- Placement
Incremental placement and light-curing each increment is strongly recommended especially in deep cavities and in Class II cavities.
- Curing
Cure the resin with a dental visible light curing unit. Hold the light tip as close to the resin as possible.

Tables: Relation between curing time and depth of cure for each dental curing unit.

Conventional halogen	(depth of cure: mm)	
	Shades	20 seconds 40 seconds
A2, A3, A3.5, B2, B3, XL	2.0	2.0
A4, HO	1.5	2.0

Fast halogen	(depth of cure: mm)	
	Shades	5 seconds 10 seconds
A2, A3, A3.5, B2, B3, XL	2.0	2.0
A4, HO	1.5	2.0

Plasma arc	(depth of cure: mm)	
	Shades	5 seconds 10 seconds
A2, A3, A3.5, B2, B3, XL	2.0	2.0
A4, HO	1.5	2.0

LED	(depth of cure: mm)	
	Shades	20 seconds 40 seconds
A2, A3, A3.5, B2, B3, XL	2.0	2.0
A4, HO	1.5	2.0

10. Finishing

Contour the restoration and adjust the occlusion using a fine diamond point. Polish with silicon rubber points or polishing discs.

[CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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CLEARFIL is a trademark of KURARAY CO., LTD. VITA is a trademark of VITA Zahnfabrik, H.Rauter GmbH & COKG Bad Sackingen, Germany.