# CLEARFIL S<sup>3</sup> BOND PLUS

[CLEARFIL TRI-S BOND PLUS]

## ENGLISH INSTRUCTIONS FOR USE

### I. INTRODUCTION

CLEARFIL TRI-S BOND PLUS is a single-component, light-cured bonding agent that allows simultaneous treatment of both dentin and enamel using one liquid.

### II. INDICATIONS

CLEARFIL TRI-S BOND PLUS is indicated for the following uses:

- Direct restorations using light-cured composite resin
- Cavity sealing as a pretreatment for indirect restorations
- Treatment of exposed root surfaces
- Intraoral repairs of fractured restorations made of porcelain, ceramic or composite resin
- Surface treatment of prosthetic restorations made of porcelain, ceramic or composite resin
- Surface treatment of dental posts
- Post cementation and core build-ups using dual-cured composite resin CLEARFIL DC CORE PLUS
- Core build-ups using light-cured composite resin

### III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers

### IV. POSSIBLE SIDE EFFECTS

The oral mucosa membrane may turn slightly white when exposed to BOND, due to the coagulation of protein. This is a temporary phenomenon that will disappear within a few days. Provide patients with precise instructions not to give any irritations with brushing to whitened parts.

### V. INCOMPATIBILITIES

- Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol can retard the curing process.
- Do not use hemostatics containing ferric compounds, since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva, due to remaining ferric ions.
- When using hemostatics containing aluminum chloride, minimize the quantity; use caution to prevent contact with the adherent surface. Failure to do so might weaken the bond strength to the tooth structure.

### VI. PRECAUTIONS

- Safety precautions**
  - This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
  - If the patient demonstrates a hypersensitivity reaction such as rash, eczema, features of inflammation, ulcer, swelling, itching, or numbness, discontinue the use of the product, remove the product and seek medical attention.
  - Use caution to prevent the product from coming in contact with the skin or getting into the eyes. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
  - If the product comes in contact with human body tissues, take the following actions:
    - <If the product gets in the eye>
      - Immediately wash the eye with copious amounts of water and consult a physician.
      - <If the product comes in contact with the skin or the oral mucosa>
        - Immediately wipe the area with a cotton pledget or a gauze moistened with alcohol; rinse with copious amounts of water.
  - Use caution to prevent the patient from accidentally swallowing the product.
  - Avoid looking directly at the dental curing light when curing the product.
  - Avoid using the same product dispensed into a well of the dispensing dish and the same applicator brush for different patients to prevent cross infection. Discard the applicator brush after use.
  - Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers.
- Handling and manipulation precautions**
  - The product should not be used for purposes other than specified in [II.INDICATIONS].
  - The use of the product is restricted to licensed dentist.
  - The product contains ethanol, a flammable substance. Do not use it near an open flame.
  - Treat appropriately in a cavity close to the pulp or accidental pulp exposure with a pulp capping agent.
  - To prevent poor performance and handling characteristic, observe the specified light curing times and other handling requirements.
  - Clean the cavity sufficiently to prevent poor bonding. If the adherent surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.
  - Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 7 minutes after dispensing. As the volatile ethanol which is contained in the BOND evaporates, the viscosity increases, thereby making it difficult to apply.
  - During the application of the BOND to the adherent surface and the duration of leaving it in place for 10 seconds after the application, move the lighting-spot out of the mouth or turn off the light to prevent the applied BOND from being exposed to the operating light because the BOND may gelate when left under such a light. In addition, after treatment for 10 seconds, dry sufficiently as soon as possible by blowing mild air.
  - BOND contains ethanol and water. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the bond does not move; otherwise, the adhesion effect will be impaired. In order to dry sufficiently, adjust the air pressure according to the shape and size of the cavity and the prosthetic appliance. Use a vacuum aspirator to prevent the bond liquid from scattering.
  - Do not use BOND in conjunction with other bonding agents. Intermixing materials may cause a change in physical properties, including a possible decrease, from the properties expected results.
  - Do not use the mixture in conjunction with other silane coupling agents except CLEARFIL PORCELAIN BOND ACTIVATOR
  - If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat again with the product.
  - If there is an adherent surface including precious metal, apply a metal adhesive primer.
  - If there is an adherent surface including uncut enamel, apply K-ETCHANT GEL to the uncut enamel, leave it in place for 10 seconds, wash with water, and then dry. After then, apply the product to an adherent surface.
  - When placing a dental post using BOND, do not use a dental core build-up material except CLEARFIL DC CORE PLUS.
  - Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
  - The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

### 3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- The product must be refrigerated (2-8°C/ 36-46°F) when not in use, and should be brought to room temperature for more than 15 minutes before using. After taking out of the refrigerator, the product must be left standing until it comes to room temperature; otherwise, the excessive amount of liquid might be depended or the liquid might ooze after use.
- The container should be lightly capped immediately after use to reduce the evaporation of the volatile solvent (ethanol contained in the product). If the liquid does not flow easily from the nozzle, the nozzle may have plugged up. Do not dispense the liquid forcibly from the plugged container.
- If the liquid may not flow easily, since it has not been used for a long term, shake the container before use.
- Keep away from extreme heat, direct sunlight or a flame.
- The product must be stored in proper places where only dental practitioners can access it.

### VII. COMPONENTS

Please see the outside of the package for contents and quantity.

#### 1) BOND

- Principal ingredients:
- 10-Methacryloyloxydecyl dihydrogen phosphate(MDP)
  - Bisphenol A diglycidylmethacrylate (Bis-GMA)
  - 2-Hydroxyethyl methacrylate (HEMA)
  - Hydrophilic aliphatic dimethacrylate
  - Hydrophobic aliphatic methacrylate
  - Colloidal silica
  - Sodium fluoride
  - dl-Camphorquinone
  - Accelerators
  - Initiators
  - Ethanol
  - Water

#### 2) Accessories

- Applicator brush (fine <silver>)
- Applicator brush (endo <black>)
- Dispensing dish
- Light blocking plate

### VIII. CLINICAL PROCEDURES

#### A. Standard procedure I (indications [1] to [3])

- Direct restorations using light-cured composite resin**
- Cavity sealing as a pretreatment for indirect restorations**
- Treatment of exposed root surfaces**

#### A-1. Moisture control

In order to produce optimal results, avoid contamination of the treatment area with saliva or blood. A rubber dam is recommended to keep the tooth clean and dry.

#### A-2. Cavity preparations

Remove any infected dentin using caries detection material (e.g. CARIES DETECTOR) as a guide and prepare the cavity in the usual manner.

#### A-3. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

#### A-4. Acid etching of uncut enamel

If there is an adherent surface including uncut enamel, apply K-ETCHANT GEL (40% phosphoric acid gel) to the uncut enamel, leave it in place for 10 seconds, wash with water, and then dry.

#### A-5. Application of BOND

- Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.
- Apply BOND to the entire cavity wall with the applicator brush. Leave it in place for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
- Dry the entire cavity wall sufficiently by blowing mild air for more than 5 seconds until the bond does not move. Use a vacuum aspirator to prevent the bond liquid from scattering.
- Light-cure BOND for 10 seconds with a dental curing unit (see table "Dental curing unit")

Table: Dental curing unit

Type	Light source	Wavelength range and light intensity
Halogen	Halogen lamp	Light intensity <sup>②</sup> of more than 300mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm
LED	Blue LED <sup>①</sup>	Light intensity <sup>③</sup> of more than 300 mW/cm <sup>2</sup> in wavelength range from 400 - 515 nm

1) Peak of emission spectrum: 450 - 480 nm

2) Evaluated according to ISO 10650-1.

3) Evaluated according to ISO 10650-2.

#### A-6. Place composite resin restorative

**A-6a. Direct restorations using light-cured composite resin**
Apply composite resins (e.g. CLEARFIL MAJESTY Esthetic, CLEARFIL MAJESTY Posterior, CLEARFIL MAJESTY Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

#### A-6b. Cavity sealing and treatment of exposed root surfaces

Apply a thin coat of composite resins (e.g. CLEARFIL MAJESTY Flow) to the tooth, and light-cure according to the manufacturer's instructions. Remove unpolymerized resin with a cotton pledget moistened with alcohol.

#### B. Standard procedure II (indications [4])

- Intraoral repairs of fractured restorations made of porcelain, ceramic or composite resin**

#### (1) When using the mixture of BOND and CLEARFIL PORCELAIN BOND ACTIVATOR for silane treatment

##### B(1)-1. Preparation of adherent surfaces

Using a diamond point, roughen the adherent surfaces, and place a bevel at the marginal area.

##### B(1)-2. Acid etching of adherent surfaces

Apply K-ETCHANT GEL to the adherent surface (including tooth structure). Leave in place for 5 seconds, wash with water, and then dry.

##### B(1)-3. Treatment of precious metal surface

If the adherent surface includes the precious metal, apply a metal adhesive primer (e.g. ALLOY PRIMER) and dry according to the manufacturer's instructions.

##### B(1)-4. Silane treatment / bonding

- Dispense one drop each of BOND and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the dispensing dish and mix them with the applicator brush.

[NOTE]

BOND and CLEARFIL PORCELAIN BOND ACTIVATOR should be mixed immediately before application.

- Apply the mixture to the entire adherent surface. Leave in place for 5 seconds. If the adherent surface includes the tooth, leave in place for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
- Dry the adherent surface sufficiently by blowing mild air for more than 5 seconds until the bond does not move. Use a vacuum aspirator to prevent the bond liquid from scattering.
- Light-cure BOND for 10 seconds with a dental curing unit (see table "Dental curing unit" in A-5.)

#### B(1)-5. Place composite resin restorative

Apply composite resin (e.g. CLEARFIL MAJESTY Esthetic, CLEARFIL MAJESTY Posterior, CLEARFIL MAJESTY Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

[NOTE]

Use an opaque resin (e.g. CLEARFIL ST OPAQUER) to mask metal color.

#### (2) When using other silane-coupling agents for silane treatment, and applying BOND on the tooth surface

##### B(2)-1. Preparation of adherent surface

Using a diamond point, roughen the adherent surface, and place a bevel at the marginal area.

##### B(2)-2. Acid etching of adherent surface

Apply K-ETCHANT GEL to the adherent surface. Leave in place for 5 seconds, wash with water, and then dry.

##### B(2)-3. Treatment of metal surface

If the adherent surface includes the metal, apply a metal adhesive primer (e.g. ALLOY PRIMER) and dry according to the manufacturer's instructions.

##### B(2)-4. Silane treatment

Apply a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER) to the adherent surface (porcelain, ceramic or composite resin) and dry according to the manufacturer's instructions.

#### B(2)-5. Bonding

- Apply BOND to the tooth structure. Leave in place for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
- Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the bond does not move. Use a vacuum aspirator to prevent the bond liquid from scattering.
- Light-cure BOND for 10 seconds with a dental curing unit (see table "Dental curing unit" in A-5.)

#### B(2)-6. Place composite resin restorative

Apply composite resins (e.g. CLEARFIL MAJESTY Esthetic, CLEARFIL MAJESTY Posterior, CLEARFIL MAJESTY Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

[NOTE]

Use an opaque resin (e.g. CLEARFIL ST OPAQUER) to mask metal color.

#### C. Standard procedure III (indications [5])

- Surface treatment of prosthetic restorations made of porcelain, ceramic or composite resin**

#### C-1. Surface preparation

As necessary, roughen the adherent surface of the prosthetic restoration by blasting with 30 to 50 µm alumina powder at air pressure of 0.1-0.4 MPa (14-58 PSI). The air pressure should be properly adjusted to suit the material and/or shape of the prosthetic restoration, using caution to prevent it from chipping. After blasting, clean it by using an ultrasonic unit for 2 minutes, followed by drying it with air stream.

#### C-2. Acid etching of adherent surface

Apply K-ETCHANT GEL to the adherent surface. Leave in place for 5 seconds, wash with water, and then dry.

#### C-3. Surface treatment

- Dispense one drop each of BOND and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the dispensing dish and mix them with the applicator brush.

[NOTE]

BOND and CLEARFIL PORCELAIN BOND ACTIVATOR should be mixed immediately before application.

- Apply the mixture to the entire adherent surface. Leave it in place for 5 seconds.

- Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the bond does not move.

#### C-4. Cementing

Cement the prosthetic restoration using resin cement according to the manufacturer's instructions.

#### D. Standard procedure IV (indications [6])

- Surface treatment of dental posts**

#### (1) When using glass fiber posts

##### D(1)-1. Acid etching of post surface

Apply K-ETCHANT GEL to the post surface. Leave it in place for 5 seconds, wash with water, and then dry.

##### D(1)-2. Surface treatment

- Dispense one drop each of BOND and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the dispensing dish and mix them with the applicator brush.
- Apply the mixture to the post surface and leave it in place for 5 seconds.
- Dry the entire post surface sufficiently by blowing mild air for more than 5 seconds until the bond does not move.

[CAUTION]

Do not blast glass fiber posts with alumina powder. Or the posts will be damaged.

#### (2) When using non-precious metal posts

##### D(2)-1. Surface preparation

Roughen the post surface by blasting with alumina powder.

##### D(2)-2. Surface treatment

- Apply BOND to the post surface with the applicator brush and leave it in place for 5 seconds.
- Dry the entire post surface sufficiently by blowing mild air for more than 5 seconds until the bond does not move.

#### E. Standard procedure V (indications [7])

- Post cementation and core build-ups using dual-cured composite resin CLEARFIL DC CORE PLUS**

#### E-1. Moisture control

In order to produce optimal results, avoid contamination of the treatment area with saliva or blood. A rubber dam is recommended to keep the tooth clean and dry.

#### E-2. Preparing root canal

Prepare and clean the root canal opening in the usual manner.

#### E-3. Post preparation

A post of proper diameter and length is fitted to the prepared root canal.

#### E-4. Treatment of post surface

##### When using glass fiber post

Treat the post surface according to the section D(1) or silane coupling agent's instructions (e.g. CLEARFIL CERAMIC PRIMER).

##### When using metal post

Treat the post surface according to the section D(2) or a metal adhesive primer's instructions (e.g. ALLOY PRIMER)

#### E-5. Application of BOND

- Apply BOND to the entire adherent surface with the applicator brush. Leave in place for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.

- Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the bond does not move. Use a vacuum aspirator to prevent the bond liquid from scattering. Remove excess BOND with a paper point. After wiping away the excess BOND, if necessary dry the adherent surface again.

[NOTE]

- To sufficiently, adjust the air pressure properly according to the shape and size of the adherent surface.

- Light-cure BOND for 10 seconds with a dental curing unit (see table "Dental curing unit" in A-5.)

#### E-6. Placing the post and core build-up

Place the post and core build-up using CLEARFIL DC CORE PLUS according to the manufacturer's instructions.

[CAUTION]

Do not use a dental core build-up material except CLEARFIL DC CORE PLUS.

#### F. Standard procedure VI (indications [8])

##### [8] Core build-ups using light-cured composite resin

Apply this method only when the adherent surface is light accessible.

#### F-1. Moisture control

In order to produce optimal results, avoid contamination of the treatment area with saliva or blood. A rubber dam is recommended to keep the tooth clean and dry.

#### F-2. Preparing root canal

Prepare and clean the root canal opening in the usual manner.

#### F-3. Post cementation

Place the post using a luting material according to the manufacturer's instructions.

#### F-4. Application of BOND

- Apply BOND to the entire adherent surface with the applicator brush. Leave in place for 10 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
- Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the bond does not move. Use a vacuum aspirator to prevent the bond liquid from scattering.
- Light-cure BOND for 10 seconds with a dental curing unit (see table "Dental curing unit" in A-5.)

#### F-5. Core build-ups

Core build-up using light-cured composite resin (e.g. CLEARFIL PHOTO CORE) according to the manufacturer's instructions.

[CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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