004 1561 R028R

CLEARFIL™ SE Protect [Bottle]

DENTAL BONDING AGENT

ENGLISH INSTRUCTIONS FOR USE

I. INTRODUCTION

CLEARFIL SE Protect is a light-cure bonding system, consisting of a self-etching primer which has an antibacterial cavity cleansing effect and a fluoride-releasing bonding agent. The primer simultaneously treats both dentin and enamel.

II. INDICATIONS

CLEARFIL SE Protect is indicated for the following uses:

- [1] Direct restorations using light-cured composite resin or compomer
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of hypersensitive teeth and/or exposed root surfaces
- [4] Intraoral repairs of fractured crowns/bridges made of porcelain, hybrid ceramics or composite resin using light-cured composite resin
- [5] Surface treatment of prosthetic appliances made of porcelain, hybrid ceramics and cured composite resin
- [6] Core build-ups using light- or dual-cured composite resin
- [7] Cavity sealing under amalgam restorations

III. CONTRAINDICATION

Patients with a history of hypersensitivity to methacrylate monomers

IV. POSSIBLE SIDE EFFECT

The oral mucosa may turn whitish when contacted by the PRIMER due to the coagulation of protein. This is a temporary phenomenon that will disappear in a few days.

V. INCOMPATIBILITIES

- [1] Do not use eugenol containing materials for pulp protection or temporary sealing, since the eugenol could retard the bonding system curing process.
- [2] Do not use hemostatics containing ferric compounds, since they may impair adhesion and the remaining ferric ion may cause discoloration in the layer or surrounding gingiva.

VI. PRECAUTIONS

1. Safety precautions

- The product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
- If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching, or numbness, discontinue the use of the product and seek medical attention.
- Avoid direct contact with the skin and/or soft tissue to prevent hypersensitivity. Wear gloves or take appropriate precautions when using the product.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
- If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye> Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin or the oral mucosa> Immediately wipe the area with a cotton pledget or a gauze moistened with alcohol; rinse with copious amounts of water.
- Use caution when using the product to prevent swallowing.
- Avoid looking directly at the dental curing light when curing the product.
- Keep the product out of reach of children.
- Avoid using the same disposable brush tip for different patients to prevent cross infection. Discard the tip after use and sterilize the brush tip handle.

2. Handling and manipulation precautions

- PRIMER and BOND will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes from dispensing.
- The emitting tip of the dental curing light should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Low intensity of light causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
- The use of the product is restricted to licensed dentists.

3. Storage precautions

- If the liquid may not flow easily since it has not been used for a long term, shake the container before use.
- The product must be used by the expiration date indicated on the package.
- The product must be refrigerated (2-8°C/36-46°F) when not in use; and should be brought to room temperature before using. After taking out of the refrigerator, the product must be left standing until it comes to room temperature; otherwise the liquid might drop more than necessary or ooze after use.
- Keep away from extreme heat or direct sunlight.
- Replace the bottle cap as soon as possible after the liquid has been dispensed from the bottle.
- The product must be stored in proper places where only dental practitioners can access it.

VII. COMPONENTS

Please see the outside of the package for contents and quantity.

1) PRIMER (self-etching primer)

- Principal ingredients:
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
 - 12-Methacryloyloxydodecylpyridinium bromide (MDPB)
 - 2-Hydroxyethyl methacrylate (HEMA)
 - Hydrophilic dimethacrylate
 - Water

2) BOND (fluoride-bonding agent)

- Principal ingredients:
- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
 - Bisphenol A diglycidylmethacrylate (Bis-GMA)
 - 2-Hydroxyethyl methacrylate (HEMA)
 - Hydrophobic dimethacrylate
 - di-Camphorquinone
 - N,N-Diethanol-p-toluidine
 - Silanated colloidal silica
 - Surface treated sodium fluoride

3) Accessories

- Disposable brush tips
- Brush tip handles
- Mixing dish
- Light blocking plate
- Outer case

[NOTE]

- * The PRIMER and BOND can easily be dispensed using the outer case container.
- * The PRIMER and BOND must be refrigerated when not in use.

VIII. CLINICAL PROCEDURES

A. Direct restorations, cavity sealing and treatment of hypersensitive teeth or exposed root surfaces using light-cured composite resin (or compomer)

A-1. Cleaning tooth structure

An adequately cleaned cavity assures maximum adhesive performance. Be sure the cavity is adequately cleaned.

A-2. Moisture control

In order to produce optimal results, avoid contamination of the treatment area from saliva or blood. A rubber dam is recommended to keep the tooth clean and dry.

A-3. Cavity preparations

Remove any infected dentin using CRIES DETECTOR as a guide and prepare the cavity in the usual manner.

[CAUTION]

The primer of this bonding system has an antibacterial property, but it is advisable to remove carious lesions thoroughly in the usual manner when preparing a cavity.

A-4. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

A-5. Acid etching uncut enamel

If there is a possibility of resin spreading over uncut enamel, apply K-ETCHANT Syringe to the enamel, let it remain for 10 seconds, wash with water, and then dry.

[CAUTION]

The use of PRIMER alone may not sufficiently condition uncut enamel. Overfilling of resin onto unetched, uncut enamel could cause marginal discoloration.

A-6. Tooth surface treatment

- Dispense the necessary amount of PRIMER into a well of the mixing dish immediately before application.
 - [CAUTION] PRIMER will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing.
- Apply PRIMER to the entire cavity wall with a disposable brush tip. Leave it in place for 20 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces for at least 20 seconds.
- After conditioning the tooth surface for 20 seconds, evaporate the volatile ingredients with a mild oil-free air stream.
 - [CAUTION] Avoid pooling of PRIMER. Do not wash after applying PRIMER. Observe the drying method and treatment time to ensure optimum adhesion. Avoid touching the treated surface. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat with PRIMER again.

A-7. Bonding

- Dispense the necessary amount of BOND into a well of the mixing dish.
 - [CAUTION] BOND will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing.
- Apply BOND to the entire surface of the cavity with a disposable brush tip.
- After applying BOND, create a uniform bond film using a gentle oil-free air flow.
 - [CAUTION] A strong air blast will scatter the bonding agent, resulting in poor adhesion.

- Light-cure the BOND for 10 seconds with a dental curing light (irradiate wave length: 400-515nm, light intensity: >300mW/cm²).

A-8. Either A-8a or A-8b

A-8a. Direct restorations using light-cured composite resin (or compomer)
Apply composite resin (e.g. CLEARFIL MAJESTY ES-2, CLEARFIL MAJESTY Posterior) or compomer into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

A-8b. Cavity sealing and treatment of hypersensitive teeth or exposed root surfaces

Apply a thin coat of flowable composite resin (e.g. CLEARFIL MAJESTY ES Flow) to the tooth, and light-cure the composite according to its' instructions for Use. Remove unpolymerized resin with a cotton pledget moistened with alcohol.

B. Intraoral repairs of fractured crowns / bridges made of porcelain, hybrid ceramics or composite resin using lightcured composite resin

B-1. Preparation of fractured surfaces

- Facing material surface
Using a diamond point, remove a layer of the fractured surface, and place a bevel at the marginal area.
- Metal surface
Roughen the metal surface with a diamond point.

B-2. Acid etching of facing material surface

Apply K-ETCHANT Syringe to the facing material surface and leave it in place for 5 seconds before washing and drying. If the adherent surface extends to uncut enamel, apply K-ETCHANT Syringe to the enamel surface and leave it in place for 10 seconds before washing and drying.

B-3. Precious metal surface

When precious metal is used, apply ALLOY PRIMER according to the manufacturer's instructions.

B-4. Silane treatment

- Dispense one drop each of PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR into a well of the mixing dish and mix them.

[NOTE]

PRIMER and CLEARFIL PORCELAIN BOND ACTIVATOR should be mixed immediately before application.

[CAUTION]

- PRIMER or the mixture will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing.
- Apply mixture to the facing material and metal surfaces with a disposable brush tip.
- Leave it in place for 5 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
- If the adherent surface extends to the tooth, leave mixture in place for 20 seconds.
- After that, evaporate the volatile ingredients with a mild oil-free air stream. The treated metal surface must be sufficiently dried; if not, adhesion might be seriously impaired.

[CAUTION]

Do not wash treated surface. If the treated surface is contaminated with saliva, wash it with water, dry, or clean with alcohol, and treat with the mixture again.

B-5. Bonding

- Dispense the necessary amount of BOND into a well of the mixing dish.
 - [CAUTION] BOND will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing.
- Apply BOND to the entire cavity wall with a disposable brush tip.
- After applying BOND, create a uniform bond film using a gentle oil-free air flow.
- Light-cure the BOND for 10 seconds with a dental curing light (irradiate wave length: 400-515nm, light intensity: >300mW/cm²).

B-6. Light-cured composite resin filling
Apply composite resin (e.g. CLEARFIL MAJESTY ES-2, CLEARFIL MAJESTY Posterior) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

[NOTE]

Use opaque resin such as CLEARFIL ST OPAQUER for metal surface to prevent metal shine through.

C. Treatment of prosthetic appliances

Silanating treatment of prosthetic appliances made of esthetic materials (porcelain, hybrid ceramics or composite resin) can be accomplished using PRIMER in conjunction with CLEARFIL PORCELAIN BOND ACTIVATOR.

C-1. Acid etching of adherent surface

Apply K-ETCHANT Syringe to the adherent surface and leave it in place for 5 seconds before washing and drying.

C-2. Silane treatment

Use the same procedure described in B-4.

C-3. Cementing

Cement the prosthetic appliance using PANAVIA cement according to the manufacturer's instructions.

D. Core build-ups using light- or dual-cured composite resin

- D-1. Cleaning tooth structure**
An adequately cleaned tooth surface assures maximum adhesive performance. Be sure the tooth surface is adequately cleaned.

D-2. Moisture control

In order to produce optimal results, avoid contamination of the treatment area from saliva or blood. A rubber dam is recommended to keep the tooth clean and dry.

D-3. Tooth preparations

Remove existing restorations and/or caries and place pin, post or matrix as needed.

D-4. Pulp protection

In case of covering the pulp exposure or near pulp exposure, a hard setting calcium hydroxide material is recommended. Do not use eugenol materials for pulp protection.

D-5. Application of primer

- Dispense the necessary amount of PRIMER into a well of the mixing dish immediately before application.

[CAUTION]

PRIMER will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing.

- Apply PRIMER to all the tooth surfaces (and the post if applicable) with a disposable brush tip. Leave it in place for 20 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces for at least 20 seconds.
- After conditioning all the tooth surfaces (and the post) for 20 seconds, evaporate the volatile ingredients with a mild oil-free air stream.

[CAUTION]

Avoid pooling of PRIMER. Do not wash after applying PRIMER. Observe the drying method and treatment time to ensure optimum adhesion. Avoid touching the treated surface. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat with PRIMER again.

D-6. Bonding

- Dispense the necessary amount of BOND into a well of the mixing dish.

[CAUTION]

BOND will set to a gel if left under an operating light or natural light (sunlight from windows). Use the light blocking plate to avoid exposing the material to an operating light or natural light, and use within 3 minutes after dispensing.

- Apply BOND to all the tooth surfaces with a disposable brush tip.
- After applying BOND, create a uniform bond film using a gentle oil-free air flow.

[CAUTION]

A strong air blast will scatter the bonding agent, resulting in poor adhesion.

D-7. Core build-ups

- Place light- or dual-cured composite resin for core-build up (e.g. CLEARFIL PHOTO CORE, CLEARFIL DC CORE PLUS), using a core instrument.
- Light-cure the core material according to the manufacturer's instructions.

E. Cavity sealing under amalgam restorations

E-1. Cleaning tooth structure

E-2. Moisture control

E-3. Cavity preparations

Prepare the cavity in the normal manner for an amalgam restoration.

E-4. Pulp protection

E-5. Acid etching uncut enamel

E-6. Tooth surface treatment

E-7. Bonding

Use the same procedures described in A-1 to A-7 for the above steps.

E-8. Amalgam filling

Mix amalgam powder and mercury, condense the mixture into the cavity, finish and polish according to the manufacture's instructions.

[CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of a dental professional.

[WARNING]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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ENGLISH INSTRUCTIONS FOR USE

I. INTRODUCTION

CLEARFIL SE Protect Unit Dose is a light-cure, self-etching bonding system, consisting of a self-etching primer which has an antibacterial cavity cleansing effect and a fluoride-releasing bonding agent. The primer simultaneously treats both dentin and enamel. The product is provided for a single-use.

II. INDICATIONS

The product is indicated for the following uses:

- [1] Direct restorations using light-cured composite resin or compomer
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of hypersensitive teeth and/or exposed root surfaces
- [4] Core build-ups using light- or dual-cured composite resin
- [5] Cavity sealing under amalgam restorations

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IV. POSSIBLE SIDE EFFECTS

The oral mucosal membrane may turn whitish when contacted by the product due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing.

V. INCOMPATIBILITIES

- [1] Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol can retard the curing process.
- [2] Do not use hemostatics containing ferric compounds since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva due to remaining ferric ions.
- [3] When using hemostatics containing aluminum chloride, minimize the quantity; use caution to prevent contact with the adherend surface. Failure to do so might weaken the bond strength to the tooth structure.

VI. PRECAUTIONS

1. Safety precautions

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components.
2. If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching or numbness, discontinue use of the product and seek medical attention.
3. Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
4. If the product comes in contact with human body tissues, take the following actions:
<If the product gets in the eye>
Immediately wash the eye with copious amounts of water and consult a physician.
<If the product comes in contact with the skin or the oral mucosa>
Immediately wipe the area with a cotton pellet or a gauze pad moistened with alcohol, and rinse with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. Avoid looking directly at the dental curing light when curing the product.
7. Single use only. Avoid using the same applicator brush for different patients to prevent cross contamination. Discard the container and applicator brush after use.
8. Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components.

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1. The product must not be used for any purposes other than specified in [II. INDICATIONS].
2. The use of this product is restricted to licensed dental professionals.
3. Use a pulp capping agent in a cavity close to the pulp or in the event of accidental pulp exposure.
4. To prevent poor performance or poor handling characteristics, observe the specified light-curing times and other handling requirements.

5. Clean the cavity sufficiently to prevent poor bonding. If the adherend surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.
6. PRIMER is applied to the adherend surface and should be left in place for 20 seconds. Move the operating light away from the mouth or turn off the light to prevent the applied PRIMER and BOND from being exposed to the operating light (the BOND may gelate when left under the light).
7. Dry sufficiently with mild air after treatment with PRIMER for 20 seconds.
8. Do not rinse after applying PRIMER.
9. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat again with the product.
10. Apply PRIMER/BOND with the applicator brush immediately after snapping off the container cap.
11. Do not use the same applicator brush for PRIMER and BOND.
12. When snapping off the container cap, do not tilt it to avoid spilling PRIMER/BOND.
13. Do not mix the product with other materials.

[Dental light-curing unit]

1. Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
2. The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface has to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

3. Storage precautions

1. The product must be used by the expiration date indicated on the package.
2. The product must be refrigerated (2-8°C/ 36-46°F) when not in use, and should be brought to room temperature for more than 15 minutes before using. After taking out of the refrigerator, the product must be left standing until it comes to room temperature.
3. The product must be kept away from extreme heat or direct sunlight.
4. The product must be stored in a proper place where only dental practitioners can access.

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- 1) PRIMER (self-etching primer) : Blue container

Principal ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- 12-Methacryloyloxydodecyl pyridinium bromide (MDPB)
- 2-Hydroxyethyl methacrylate
- Hydrophilic dimethacrylate
- Water

- 2) BOND (fluoride-bonding agent) : Black container

Principal ingredients:

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
- Bisphenol A diglycidylmethacrylate
- 2-Hydroxyethyl methacrylate
- Hydrophobic dimethacrylate
- di-Camphorquinone
- N,N-Diethanol-p-toluidine
- Silanated colloidal silica
- Surface treated sodium fluoride

- 3) Accessories

Applicator brush (fine <silver>)

VIII. CLINICAL PROCEDURES

A. Standard procedure I

- [1] Direct restorations using light-cured composite resin or compomer
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of hypersensitive and/or exposed root surfaces

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A-3. Cavity preparations

Remove any infected dentin using CRIES DETECTOR as a guide and prepare the cavity in the usual manner. [CAUTION]

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A-4. Pulp protection

Any actual or near pulp exposure can be covered with a hard setting calcium hydroxide material. There is no need for a cement liner or base. Do not use eugenol materials for pulp protection.

A-5. Acid etching uncut enamel

If there is an adherend surface including uncut enamel, apply a phosphoric acid-etching gel (e.g. K-ETCHANT Syringe) to the uncut enamel, leave it in place for 10 seconds, then rinse and dry.

A-6. Application of PRIMER

1. Snap off the container cap and apply PRIMER to the entire cavity wall with the applicator brush, and then leave it in place for 20 seconds. Use caution not to allow saliva or exudate to contact the treated surfaces.
2. Dry the entire cavity wall sufficiently with mild air for more than 5 seconds until the PRIMER does not move. Use a vacuum aspirator to prevent the PRIMER from scattering. [CAUTION]
Do not rinse after applying PRIMER.

A-7. Application of BOND

1. Snap off the container cap and apply BOND to the entire cavity wall with an applicator brush. Use caution not to allow saliva or exudate to contact the treated surfaces.
2. Make a uniform bond film using a gentle air flow.
3. Light-cure BOND for 10 seconds with a dental curing unit (irradiate wave length: 400-515nm, light intensity: >300mW/cm²).

A-8. Placement of composite resin restorative or Treatment of hypersensitive teeth

A-8a. Direct restorations using light-cured composite resin or compomer

Apply composite resins (e.g. CLEARFIL MAJESTY ES-2, CLEARFIL MAJESTY ES Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

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Apply a thin coat of composite resins (e.g. CLEARFIL MAJESTY ES Flow) onto the tooth, and light-cure according to the manufacturer's instructions. Remove unpolymerized resin with a cotton pellet or a gauze moistened with alcohol.

B. Standard procedure II

- [4] Core build-ups using light- or dual-cured composite resin

B-1. Cleaning tooth structure

An adequately cleaned tooth surface assures maximum adhesive performance. Be sure the tooth surface is adequately cleaned.

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B-5. Application of PRIMER and BOND

Use the same procedure described in "A-6" and "A-7".

B-6. Core build-up

Perform a core build-up using a material (e.g. CLEARFIL PHOTO CORE, CLEARFIL DC CORE PLUS) according to the manufacturer's instructions.

[CAUTION]

When using a dual-cure core material, the first increment should be light-cured within the depth of cure of the material.

C. Standard Procedure III

- [5] Cavity sealing under amalgam restorations

C-1. Cleaning tooth structure, moisture control

Use the same procedure described in "A-1" and "A-2".

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C-3. Pulp protection

Use the same procedure described in "A-4".

C-4. Acid etching uncut enamel

Use the same procedure described in "A-5".

C-5. Application of PRIMER and BOND

Use the same procedure described in "A-6" and "A-7".

C-6. Amalgam filling

Mix amalgam powder and mercury, condense the mixture into the cavity, finish and polish according to the manufacturer's instructions.

[CAUTION]

Federal (U.S.A.) law restricts this device to sale by or on the order of dental professionals.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

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"CLEARFIL", "CLEARFIL SE Protect", "CLEARFIL MAJESTY", "PANAVIA" and "CLEARFIL DC CORE PLUS" are trademarks of KURARAY CO., LTD.
"CLEARFIL", "CLEARFIL MAJESTY", "CLEARFIL DC CORE PLUS", "CLEARFIL PHOTO CORE" and "CLEARFIL SE Protect" are trademarks of KURARAY CO., LTD.
La legislación federal (U.S.A.) limita este dispositivo a la venta por o bajo pedido de un profesionales dentales.