

Safety Data Sheet
according to HPR, Schedule 1

Printing date 09/25/2023

Version 4

Reviewed on 09/25/2023

1 Identification

- **Product identifier**
- **Trade name:** CLEARFIL SE PROTECT ; BOND
- **Article number:** CA092-B
- **Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the mixture** Dental bonding material
- **Details of the supplier of the safety data sheet**
- **Supplier**
Kuraray America, Inc.
32 Old Slip, 7th Floor, New York,
NY 10005 U.S.A.
Tel: 800-879-1676
Website: www.kuraraydental.com
- **Further information obtainable from:** Supplier
- **Emergency telephone number:**
CHEMTREC; day and night, within USA and Canada
1-800-424-9300 CCN706984 or +1 703-527-3887 (collect calls accepted)

2 Hazard identification

- **Classification of the substance or mixture**
Skin Irritation - Category 2 H315 Causes skin irritation.
Eye Irritation - Category 2A H319 Causes serious eye irritation.
Skin Sensitizer - Category 1 H317 May cause an allergic skin reaction.
Reproductive Toxicity - Category 2 H361 Suspected of damaging fertility or the unborn child.

- **Label elements**
- **GHS label elements**
The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms**



GHS07 GHS08

- **Signal word** Warning
- **Hazard-determining components of labeling:**
2-hydroxyethyl methacrylate
diphenyl(2,4,6- trimethylbenzoyl)phosphine oxide
- **Hazard statements**
H315 Causes skin irritation.
H319 Causes serious eye irritation.
H317 May cause an allergic skin reaction.
H361 Suspected of damaging fertility or the unborn child.
- **Precautionary statements**
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313 IF exposed or concerned: Get medical advice/attention.
P321 Specific treatment (see on this label).
P362+P364 Take off contaminated clothing and wash it before reuse.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P337+P313 If eye irritation persists: Get medical advice/attention.

(Contd. on page 2)

Safety Data Sheet

according to HPR, Schedule 1

Printing date 09/25/2023

Version 4

Reviewed on 09/25/2023

Trade name: CLEARFIL SE PROTECT ; BOND

P405 Store locked up.
 P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

(Contd. of page 1)

3 Composition/Information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of substances listed below with nonhazardous additions.

· **Dangerous components:**

1565-94-2	Bisphenol A diglycidylmethacrylate Skin Irritation - Category 2, H315; Eye Irritation - Category 2A, H319	25-45% *
868-77-9	2-hydroxyethyl methacrylate Skin Irritation - Category 2, H315; Eye Irritation - Category 2A, H319; Skin Sensitizer - Category 1, H317	20-40% *
75980-60-8	diphenyl(2,4,6- trimethylbenzoyl)phosphine oxide Reproductive Toxicity - Category 2, H361	<3%
7681-49-4	sodium fluoride Acute Toxicity (Oral) - Category 3, H301; Skin Irritation - Category 2, H315; Eye Irritation - Category 2A, H319	<1%

* Actual concentration ranges are withheld as a trade secret.

· **Other ingredients:**

Hydrophobic aliphatic dimethacrylate
 Amorphous silica
 10-Methacryloyloxydecyl dihydrogen phosphate
 dl-Camphorquinone
 Accelerator

4 First-aid measures

- **Description of first aid measures**
- **After inhalation:**
 Supply fresh air and to be sure call for a doctor.
 In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:**
 Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:**
 Do not induce vomiting; immediately call for medical help.
 If symptoms persist consult doctor.
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
 No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
 CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:**
 Wear self-contained respiratory protective device.

(Contd. on page 3)

Safety Data Sheet
according to HPR, Schedule 1

Printing date 09/25/2023

Version 4

Reviewed on 09/25/2023

Trade name: CLEARFIL SE PROTECT ; BOND

Wear fully protective suit. (Contd. of page 2)

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** See Section 8.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and storage

- **Precautions for safe handling**
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store in a refrigerator (2-8°C).
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**
Keep receptacle tightly sealed.
Store in cool, dry conditions in well sealed receptacles.
Protect from heat and direct sunlight.
- **Specific end use(s)** No further relevant information available.

8 Exposure controls/ Personal protection

- **Control parameters**

Components with limit values that require monitoring at the workplace:	
7681-49-4 sodium fluoride	
EL	TWA: 2.5 mg/m ³ as F

- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Additional information about design of technical systems:** No further data; see section 7.
- **Personal protective equipment:**
- **General protective and hygienic measures:**
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.
- **Breathing equipment:** Not required.
- **Protection of hands:**



Protective gloves

Neoprene gloves

(Contd. on page 4)

Trade name: CLEARFIL SE PROTECT ; BOND

(Contd. of page 3)

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

- **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

- **Penetration time of glove material**

The exact break trough time has to be found out by the manufacturer of the protective gloves and has to be observed.

- **Eye protection:**



Tightly sealed goggles

9 Physical and chemical properties

- **Information on basic physical and chemical properties**

· Physical state	Liquid
· Color:	Transparent Yellow
· Odor:	Odorless
· Melting point/Melting range:	-12 °C (2-hydroxyethyl methacrylate)
· Boiling point/Boiling range:	205 °C (2-hydroxyethyl methacrylate)
· Flammability (solid, gaseous):	Not determined.
· Explosion limits:	
· Lower:	Not determined
· Upper:	Not determined
· Flash point:	97 °C (2-hydroxyethyl methacrylate)
· Auto ignition temperature:	Not determined
· Decomposition temperature:	Not determined
· pH-value:	Not determined
· Viscosity:	Not determined
· Kinematic viscosity	Not determined.
· Solubility in / Miscibility with	
· Water:	Not miscible or difficult to mix.
· Partition coefficient (n-octanol/water):	Not determined.
· Vapor pressure:	Not determined
· Density:	1.1 g/cm ³
· Relative density	Not determined
· Relative vapor density	Not determined

- **Other information**

- **Important information on protection of health and environment, and on safety.**

- **Danger of explosion:** Product does not present an explosion hazard.

10 Stability and reactivity

- **Reactivity** No further relevant information available.

(Contd. on page 5)

Safety Data Sheet
according to HPR, Schedule 1

Printing date 09/25/2023

Version 4

Reviewed on 09/25/2023

Trade name: CLEARFIL SE PROTECT ; BOND

(Contd. of page 4)

- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
No decomposition if used and stored according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** Carbon monoxide and carbon dioxide

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:** Based on available data, the classification criteria are not met.

LD/LC50 values that are relevant for classification:		
868-77-9 2-hydroxyethyl methacrylate		
Oral	LD50	4680 mg/kg (guinea pig) 3275 mg/kg (mouse) 5050 mg/kg (rat)
7681-49-4 sodium fluoride		
Oral	LD50	44 mg/kg (mouse) 52 mg/kg (rat)

- **on the skin:**
Causes skin irritation.
- **on the eye:**
Causes serious eye irritation.
- **Sensitization:**
May cause an allergic skin reaction.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity**
Suspected of damaging fertility or the unborn child.
- **Specific target organ toxicity - single exposure**
Based on available data, the classification criteria are not met.
- **Specific target organ toxicity - repeated exposure**
Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **Subacute to chronic toxicity:**

· Carcinogenic categories
· IARC (International Agency for Research on Cancer)
None of the ingredients is listed.

· NTP (National Toxicology Program)
None of the ingredients is listed.

12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.

(Contd. on page 6)

Safety Data Sheet
according to HPR, Schedule 1

Printing date 09/25/2023

Version 4

Reviewed on 09/25/2023

Trade name: CLEARFIL SE PROTECT ; BOND

(Contd. of page 5)

- **vPvB:** Not applicable.
- **Other adverse effects**
- **Additional ecological information:**
- **General notes:**
Water hazard class 1 (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

13 Disposal considerations

- **Recommendation:**
Must not be disposed of together with household garbage. Do not allow product to reach sewage system.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

14 Transport information

· UN-Number	Void
· DOT/TDG, IMDG, IATA	Void
· UN proper shipping name	Void
· DOT/TDG, IMDG, IATA	Void
· Transport hazard class(es)	
· DOT, IMDG, IATA	
· Class	Void
· Packing group	
· DOT/TDG, IMDG, IATA	Void
· Environmental hazards:	
· Marine pollutant:	No
· Special precautions for user	Not applicable.
· Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code	Not applicable.
· UN "Model Regulation":	Void

15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**
No further relevant information available.

· **Canadian substance listings:**

· Canadian Domestic Substances List (DSL)	
1565-94-2	Bisphenol A diglycidylmethacrylate
868-77-9	2-hydroxyethyl methacrylate
7681-49-4	sodium fluoride
· Canadian Non-Domestic Substances List (NDSL)	
None of the ingredients is listed.	
· Canadian Ingredient Disclosure list (limit 0.1%)	
None of the ingredients is listed.	

(Contd. on page 7)

Printing date 09/25/2023

Version 4

Reviewed on 09/25/2023

Trade name: CLEARFIL SE PROTECT ; BOND

(Contd. of page 6)

· **Canadian Ingredient Disclosure list (limit 1%)**

868-77-9	2-hydroxyethyl methacrylate
7681-49-4	sodium fluoride

16 Other information

· **Relevant phrases**

- H301 Toxic if swallowed.
- H315 Causes skin irritation.
- H317 May cause an allergic skin reaction.
- H319 Causes serious eye irritation.
- H361 Suspected of damaging fertility or the unborn child.

· **Abbreviations and acronyms:**

- DOT: US Department of Transportation
- ELINCS: European List of Notified Chemical Substances

· *** Data compared to the previous version altered.**

CA