

		Page Number	Classification	Application Time	Curing Type	Curing time (light-cure)	Direct Restoration	Indirect Restoration	Combination with Self-cure Composite	Desensitization	Film Thickness
Universal											
CLEARFIL™ Universal Bond Quick		14	1-step system	3 sec. (with no wait)	Light	10 sec.	✓	✓**		✓	5-10 microns
Self-etch											
CLEARFIL™ SE BOND 2		16	2-step system	20 sec.	Light	10 sec.	✓	✓**		✓	50 microns
CLEARFIL™ SE BOND		18	2-step system	20 sec.	Light	10 sec.	✓	✓**		✓	50 microns
CLEARFIL™ SE Protect		20	2-step system	20 sec.	Light	10 sec.	✓	✓**		✓	50 microns
CLEARFIL™ LINER BOND 2V		22	2-step system	30 sec.	Light / Dual*	10 sec. (optional)	✓	✓	✓	✓	50 microns
Total-etch											
CLEARFIL™ PHOTO BOND		23	2-step system (including etching)	120 sec.	Dual	10 sec.	✓	✓	✓		1 micron
CLEARFIL™ NEW BOND		24	2-step system (including etching)	60 sec.	Self	N/A	✓	✓	✓		1 micron

* Dual-cure is achieved when Bond A is mixed with Bond B.

** with use of CLEARFIL DC Activator.