










	Indications					Adhesion to					Curing systems		
	Crowns & bridges	Inlays/onlays	Veneers	Maryland bridges	Posts	Metal oxide ceramics	Silica-based ceramics	Precious metal	Non-precious metal	Tooth Primer	Curing type	Fluoride release	Tack Cure Clean-up
 <p>PANAVIA™ SA Cement Universal</p>  	✓	✓		✓	✓	✓	✓	✓ 2)	✓	silane built into paste	dual-cure	✓	✓
 <p>PANAVIA™ V5</p> 	✓	✓	✓ (best for esthetics)	✓	✓ 1)	✓ 1)	✓ 1)	✓ 2)	✓ 1)	✓	dual-cure	✓	✓
 <p>PANAVIA™ F 2.0</p> 	✓	✓	✓	✓	✓	✓	✓ 3)	✓ 2)	✓	✓	dual-cure	✓	
 <p>PANAVIA™ 21</p> 	✓	✓		✓	✓	✓	✓ 3)	✓ 2)	✓		self-cure		

1) Must be used with CLEARFIL™ CERAMIC PRIMER PLUS 2) Must be used with ALLOY PRIMER
3) Must be used with CLEARFIL™ CERAMIC PRIMER PLUS or CLEARFIL™ PORCELAIN BOND ACTIVATOR