

kuraray

DENTAL LIGHT-CURED RESTORATIVE COMPOSITE

CLEARFIL MAJESTY™ ES-2

A Universal Nano-Hybrid Composite

Intuitive Shade Matching—with a Single Shade

- **New: The ONLY composite with VITA Approved Shades** which match the VITA Classical Shade Guide
- **Unsurpassed natural shade blending** with surrounding tooth structure
- **Single shade placement technique**, no complex multi-shade layering due to Kuraray's pioneering LD (Light-Diffusion) Technology
- **Easy polishing** with durable gloss
- **Long working time** for easy sculpting with no slumping
- **Universal** for anterior and posterior use
- **Virtually no post-cure color change**



KURARAY AMERICA, INC.

33 Maiden Ln, Suite 600-D, New York, NY 10038 | 800.879.1676 | kuraraydental.com | ©2014 Kuraray America, Inc.



Single Shade, intuitive color matching with Kuraray's New LD Technology - a light diffusion technology that is like having millions of micro-sized prisms in the composite.

The intuitive shade matching feature of CLEARFIL™ MAJESTY ES-2 comes from Kuraray's proprietary new LD Technology, a light diffusion technology that is like millions of micro-sized prisms all refracting and transmitting color and light from the surrounding dentin and enamel through the composite. Precise color blending that begins at placement and continues from one visit to the next.

Shade Matching Like Never Before With Kuraray's LD Technology!


This series of images shows how Kuraray's LD Technology works.

Natural Tooth	Restoration with Conventional Body Shade only	Restoration with Conventional Enamel Shade only	Restoration with Conventional Layer Technique	Restoration with Single-Shade CLEARFIL MAJESTY™ ES-2
Light enters through enamel, passes in straight lines, but diffuses through dentin.	Light bounces back resulting in unnatural looking tooth.	Light passes through completely but the intraoral darkness causes a grayish, lifeless color.	Layering using both enamel and dentin shades can result in a more natural restoration.	Light enters the surface and is diffused as it penetrates deeper, just like dentin. CLEARFIL MAJESTY™ ES-2 replicates natural tooth structure.

CLEARFIL MAJESTY™ ES-2 restorations are virtually invisible due to Kuraray's LD Technology.



Photos courtesy of Dr. Hirofumi Tashiro



Warranty: All data, other than from quoted articles herein, are based on actual measurements performed by Kuraray Noritake Dental Inc. Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Neither Kuraray Noritake Dental Inc. nor the distributor accepts liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user should determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therein. CLEARFIL MAJESTY™ ES-2 and Kuraray Noritake Dental Inc. products are available only through Kuraray Authorized dealers.

ORDERING INFORMATION

CLEARFIL MAJESTY™ ES-2

Classic

Shade	Syringe 2.0ml/3.6g	PLT 0.25gx20	Shade	Syringe 2.0ml/3.6g	PLT 0.25gx20
A1	#3020KA	#3110KA	B4	#3029KA	#3119KA
A2	#3021KA	#3111KA	C1	#3030KA	#3120KA
A3	#3022KA	#3112KA	C2	#3031KA	#3121KA
A3.5	#3023KA	#3113KA	C3	#3032KA	#3122KA
A4	#3024KA	#3114KA	C4	#3033KA	#3123KA
A6	#3025KA	#3115KA	D2	#3034KA	#3124KA
B1	#3026KA	#3116KA	D3	#3035KA	#3125KA
B2	#3027KA	#3117KA	XW	#3036KA	#3126KA
B3	#3028KA	#3118KA	W	#3037KA	#3127KA

ACCESSORY

#3000KA CLEARFIL MAJESTY ES-2 SHADE GUIDE