

Recommended
 Adequate
 OK, but not recommended

	Page Number	Class I	Class II	Class III	Class IV	Class V	Flexural Strength	Filler	Volumetric Shrinkage	Radiopaque	Shades
	44	 ■	 ■	 ■	 ■	 ■	118MPa	78wt%	1.9%	✓	18 VITA Approved Classic Shades 5 Premium Shades (Enamel & Dentin) 4 Translucent Shades
	48	■	■	■	■*	■	151MPa	75wt%	2.9%	✓	10 shades
	49	■	■	■	■*	■	145MPa	81wt%	2.9%	✓	12 shades
	50	■	■	■	■	■	177MPa	92wt%	1.5%	✓	6 shades
	51	■	■	■	■	■	164MPa	85.5wt%	1.9%	✓	14 shades for Syringes, 8 shades for PLTs

*for liner use only