




# KATANA™ Zirconia Block

**Compatible With Block Scanner  
for CEREC Primemill™**

## Automatically Scans The Zirconia Block

In order to support the block scanner system introduced in the new CEREC Primemill™ system sold by Dentsply Sirona, the traditional two-dimensional barcode has been replaced by a data matrix code.

### Details



Existing CEREC users who own MC series units (prior to Primemill) can utilize the new data matrix code to initiate the touch process.

There will be a seven alphanumeric code displayed on each block. This will resolve previous “Quality Label Error” codes which required entering a space after only six characters.



The printed side of the bar code has been changed to reflect the new data matrix code.



The printed side of the KATANA™ logo remains unchanged.



**Kuraray Noritake Dental Inc.**

300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan

Website [www.kuraraynoritake.com](http://www.kuraraynoritake.com)

- Before using this product, be sure to read the Instructions for Use supplied with the product.
- The specifications and appearance of the product are subject to change without notice.
- Printed color can be slightly different from actual color.

“KATANA” is trademark of NORITAKE CO.,LIMITED.

“CEREC Primemill” is trademark of Dentsply Sirona Inc.