

Safety Data Sheet
according to HPR, Schedule 1

Printing date 10/26/2023

Version 4

Reviewed on 10/26/2023

1 Identification

- **Product identifier**
- **Trade name:** CLEARFIL PORCELAIN BOND ACTIVATOR
- **Article number:** 028
- **Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the mixture** Dental materials
- **Details of the supplier of the safety data sheet**
- **Supplier**
Kuraray America, Inc.
32 Old Slip, 7th Floor, New York,
NY 10005 U.S.A.
Tel: 800-879-1676
Website: www.kuraraydental.com
- **Further information obtainable from:** Supplier
- **Emergency telephone number:**
CHEMTREC; day and night, within USA and Canada
1-800-424-9300 CCN706984 or +1 703-527-3887 (collect calls accepted)

2 Hazard identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonized System (GHS).
- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void

3 Composition/Information on ingredients

- **Chemical characterization:** Mixtures
- **Description:** Mixture of substances listed below with nonhazardous additions.
- **Dangerous components:** Void
- **Other ingredients:**
Hydrophobic aromatic dimethacrylate
3-Trimethoxysilylpropyl methacrylate

4 First-aid measures

- **Description of first aid measures**
- **After inhalation:**
Supply fresh air and to be sure call for a doctor.
In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:**
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:**
Do not induce vomiting; immediately call for medical help.
If symptoms persist consult doctor.
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.

(Contd. on page 2)

CA

Safety Data Sheet
according to HPR, Schedule 1

Printing date 10/26/2023

Version 4

Reviewed on 10/26/2023

Trade name: CLEARFIL PORCELAIN BOND ACTIVATOR

(Contd. of page 1)

- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:**
Wear self-contained respiratory protective device.
Wear fully protective suit.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** See Section 8.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Ensure adequate ventilation.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and storage

- **Precautions for safe handling**
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store in a refrigerator (2-8°C).
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**
Keep receptacle tightly sealed.
Store in cool, dry conditions in well sealed receptacles.
Protect from heat and direct sunlight.
- **Specific end use(s)** No further relevant information available.

*** 8 Exposure controls/ Personal protection**

- **Control parameters**
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Additional information about design of technical systems:** No further data; see section 7.

(Contd. on page 3)

CA

Safety Data Sheet
according to HPR, Schedule 1

Printing date 10/26/2023

Version 4

Reviewed on 10/26/2023

Trade name: CLEARFIL PORCELAIN BOND ACTIVATOR

(Contd. of page 2)

- **Personal protective equipment:**
- **General protective and hygienic measures:**
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.
- **Breathing equipment:**
In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Safety glasses

*

9 Physical and chemical properties

- **Information on basic physical and chemical properties**
- **Physical state** Liquid
- **Color:** Transparent
- **Odor:** Slightly methacrylate odour
- **Melting point/Melting range:** Not determined
- **Boiling point/Boiling range:** 190 °C (3-trimethoxysilylpropyl methacrylate)
- **Flammability (solid, gaseous):** Not determined.
- **Explosion limits:**
- **Lower:** Not determined
- **Upper:** Not determined
- **Flash point:** 92 °C (3-trimethoxysilylpropyl methacrylate)
- **Auto ignition temperature:** Not applicable
- **Decomposition temperature:** Not determined
- **pH-value:** Not determined
- **Viscosity:** Not determined
- **Kinematic viscosity** Not determined.
- **Solubility in / Miscibility with**
- **Water:** Not miscible or difficult to mix.
- **Partition coefficient (n-octanol/water):** Not determined.
- **Vapor pressure:**
- **Density:** 1.1 g/cm³
1.8 g/cm³
- **Relative density** Not determined
- **Relative vapor density** Not determined

(Contd. on page 4)

CA

Safety Data Sheet
according to HPR, Schedule 1

Printing date 10/26/2023

Version 4

Reviewed on 10/26/2023

Trade name: CLEARFIL PORCELAIN BOND ACTIVATOR

(Contd. of page 3)

- **Other information**
- **Important information on protection of health and environment, and on safety.**
- **Danger of explosion:** Product does not present an explosion hazard.

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
No decomposition if used and stored according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** Carbon monoxide and carbon dioxide

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:** Based on available data, the classification criteria are not met.
- **LD/LC50 values that are relevant for classification:**

2530-85-0 3-trimethoxysilylpropyl methacrylate		
Oral	LD50	22,600 mg/kg (rat)
Dermal	LD50	500 mg/kg (rabbit)
- **on the skin:** Based on available data, the classification criteria are not met.
- **on the eye:** Based on available data, the classification criteria are not met.
- **Sensitization:** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **Specific target organ toxicity - single exposure**
Based on available data, the classification criteria are not met.
- **Specific target organ toxicity - repeated exposure**
Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **Subacute to chronic toxicity:**
- **Carcinogenic categories**

- **IARC (International Agency for Research on Cancer)**

None of the ingredients is listed.

- **NTP (National Toxicology Program)**

None of the ingredients is listed.

12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.

(Contd. on page 5)

CA

Safety Data Sheet
according to HPR, Schedule 1

Printing date 10/26/2023

Version 4

Reviewed on 10/26/2023

Trade name: CLEARFIL PORCELAIN BOND ACTIVATOR

(Contd. of page 4)

- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
Water hazard class 1 (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

13 Disposal considerations

- **Recommendation:**
Must not be disposed of together with household garbage. Do not allow product to reach sewage system.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

14 Transport information

- | | |
|---|-----------------|
| · UN-Number | |
| · DOT/TDG | Void |
| · IMDG, IATA | Void |
| · UN proper shipping name | |
| · DOT/TDG, IMDG, IATA | Void |
| · Transport hazard class(es) | |
| · DOT, IMDG, IATA | |
| · Class | Void |
| · Packing group | |
| · DOT/TDG | Void |
| · IMDG, IATA | Void |
| · Environmental hazards: | |
| · Marine pollutant: | No |
| · Special precautions for user | Not applicable. |
| · Transport in bulk according to Annex II of MARPOL/73/78 and the IBC Code | Not applicable. |
| · UN "Model Regulation": | Void |

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Canadian substance listings:**

· Canadian Domestic Substances List (DSL)
2530-85-0 3-trimethoxysilylpropyl methacrylate

· Canadian Non-Domestic Substances List (NDSL)
None of the ingredients is listed.

· Canadian Ingredient Disclosure list (limit 0.1%)
None of the ingredients is listed.

(Contd. on page 6)

CA

Safety Data Sheet

according to HPR, Schedule 1

Printing date 10/26/2023

Version 4

Reviewed on 10/26/2023

Trade name: CLEARFIL PORCELAIN BOND ACTIVATOR

(Contd. of page 5)

· Canadian Ingredient Disclosure list (limit 1%)

2530-85-0	3-trimethoxysilylpropyl methacrylate
-----------	--------------------------------------

*

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Abbreviations and acronyms:

DOT: US Department of Transportation

ELINCS: European List of Notified Chemical Substances

· * Data compared to the previous version altered.

CA