

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 10/27/2023

Version 6

Reviewed on 10/27/2023

**1 Identification**

- **Product identifier**
- **Trade name:** PANAVIA 21 ; OXYGUARD II
- **Article number:** US050-O
- **Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the mixture** Dental materials
- **Details of the supplier of the safety data sheet**
- **Manufacturer:**  
Kuraray Noritake Dental Inc.  
1621 Sakazu, Kurashiki,  
Okayama 710-0801,  
Japan
- **Supplier**  
Kuraray America, Inc.  
32 Old Slip, 7th Floor, New York,  
NY 10005 U.S.A.  
Tel: 800-879-1676  
Website: www.kuraraydental.com
- **Emergency telephone number:**  
CHEMTREC; day and night, within USA and Canada  
1-800-424-9300 CCN706984 or +1 703-527-3887 (collect calls accepted)

**2 Hazard(s) identification**

- **Classification of the substance or mixture**  
The product is not classified, according to the Globally Harmonized System (GHS).
- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Other hazards** Avoid breathing sanding dust / swarf.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

**3 Composition/information on ingredients**

- **Chemical characterization:** Mixtures
- **Description:** Mixture of substances listed below with nonhazardous additions.
- **Dangerous components:** Void
- **Other ingredients:**  
Glycerol  
Polyethyleneglycol  
Catalyst  
Accelerator  
Dye

**4 First-aid measures**

- **Description of first aid measures**
- **General information:**  
Avoid contact with eyes.

(Contd. on page 2)

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 10/27/2023

Version 6

Reviewed on 10/27/2023

**Trade name: PANA VIA 21 ; OXYGUARD II**

(Contd. of page 1)

- Avoid contact with skin.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If symptoms persist consult doctor.
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

**5 Fire-fighting measures**

- **Extinguishing media**
- **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

**6 Accidental release measures**

- **Personal precautions, protective equipment and emergency procedures** See Section 8.
- **Environmental precautions:**  
Dilute with plenty of water.  
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**  
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**  
No dangerous substances are released.  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

**7 Handling and storage**

- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store in a refrigerator (2-8°C).
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**  
Store in cool, dry conditions in well sealed receptacles.  
Keep receptacle tightly sealed.  
Protect from heat and direct sunlight.
- **Specific end use(s)** No further relevant information available.

**8 Exposure controls/personal protection**

- **Control parameters**

- **Components with limit values that require monitoring at the workplace:**

**56-81-5 Glycerol**

PEL	Long-term value: 15* 5** mg/m <sup>3</sup> mist; *total dust **respirable fraction
-----	---

(Contd. on page 3)

US

# Safety Data Sheet

acc. to OSHA HCS

Printing date 10/27/2023

Version 6

Reviewed on 10/27/2023

Trade name: PANA VIA 21 ; OXYGUARD II

(Contd. of page 2)

TLV	TLV withdrawn-insufficient data human occup. exp.
<b>25322-68-3 Polyethylene glycol</b>	
WEEL	Long-term value: 10 mg/m <sup>3</sup> (H); MW>200

· **Additional information:** The lists that were valid during the creation were used as basis.

· **Exposure controls**

· **Additional information about design of technical systems:** No further data; see section 7.

· **Personal protective equipment:**

· **General protective and hygienic measures:**

The usual precautionary measures for handling chemicals should be followed.

· **Breathing equipment:** Not required.

· **Protection of hands:**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:** Safety glasses

## 9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· <b>Physical state</b>	Liquid
· <b>Color:</b>	Green
· <b>Odor:</b>	Odorless
· <b>Odor threshold:</b>	Not determined.
· <b>Melting point/Melting range:</b>	18.2 °C (glycerol)
· <b>Boiling point/Boiling range:</b>	290.0 °C (glycerol)
· <b>Flammability (solid, gaseous):</b>	Not applicable.
· <b>Explosion limits:</b>	
· <b>Lower:</b>	0.9 Vol % (glycerol)
· <b>Upper:</b>	Not determined.
· <b>Flash point:</b>	160 °C (glycerol)
· <b>Auto ignition temperature:</b>	400.0 °C (glycerol)
· <b>Decomposition temperature:</b>	Not determined
· <b>pH-value:</b>	5.8 - 6.2
· <b>Viscosity:</b>	Not determined
· <b>Kinematic viscosity</b>	Not determined
· <b>Dynamic:</b>	Not determined
· <b>Solubility in / Miscibility with</b>	
· <b>Water:</b>	Fully miscible.
· <b>Partition coefficient (n-octanol/water):</b>	Not determined.
· <b>Vapor pressure:</b>	0.1 hPa (20°C, glycerol)
· <b>Density:</b>	Not determined.
· <b>Relative density</b>	Not determined

(Contd. on page 4)

US

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 10/27/2023

Version 6

Reviewed on 10/27/2023

**Trade name: PANA VIA 21 ; OXYGUARD II**

(Contd. of page 3)

- |  |   |
|--|---|
| · <b>Relative vapor density</b>  | Not determined.                               |
| · <b>Other information</b>   |   |
| · <b>Important information on protection of health and environment, and on safety.</b> |   |
| · <b>Danger of explosion:</b>  | Product does not present an explosion hazard. |
| · <b>Change in condition</b>   |   |
| · <b>Evaporation rate</b>  | Not determined.                               |
|  | Not determined                                |

**10 Stability and reactivity**

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:**  
No decomposition if used and stored according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** Carbon monoxide and carbon dioxide

**11 Toxicological information**

- **Information on toxicological effects**
- **Acute toxicity:** Based on available data, the classification criteria are not met.
- **on the skin:** Based on available data, the classification criteria are not met.
- **on the eye:** Based on available data, the classification criteria are not met.
- **Sensitization:** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **Specific target organ toxicity - single exposure**  
Based on available data, the classification criteria are not met.
- **Specific target organ toxicity - repeated exposure**  
Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **Subacute to chronic toxicity:**

- **Carcinogenic categories**

**· IARC (International Agency for Research on Cancer)**

None of the ingredients is listed.

**· NTP (National Toxicology Program)**

None of the ingredients is listed.

**12 Ecological information**

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.

(Contd. on page 5)

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 10/27/2023

Version 6

Reviewed on 10/27/2023

**Trade name: PANA VIA 21 ; OXYGUARD II**

(Contd. of page 4)

- **vPvB:** Not applicable.
- **Other adverse effects**
- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (Self-assessment): slightly hazardous for water  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

**13 Disposal considerations**

- **Recommendation:** Smaller quantities can be disposed of with household waste.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

**14 Transport information**

- |  |                 |
|--|-----------------|
| · <b>UN-Number</b>   |                 |
| · <b>DOT, IMDG, IATA</b>   | Void            |
| · <b>UN proper shipping name</b>   |                 |
| · <b>DOT, IMDG, IATA</b>   | Void            |
| · <b>Transport hazard class(es)</b>  |                 |
| · <b>DOT, IMDG, IATA</b>   |                 |
| · <b>Class</b>   | Void            |
| · <b>Packing group</b>   |                 |
| · <b>DOT, IMDG, IATA</b>   | Void            |
| · <b>Environmental hazards:</b>  |                 |
| · <b>Marine pollutant:</b>   | No              |
| · <b>Special precautions for user</b>  | Not applicable. |
| · <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b> | Not applicable. |
| · <b>UN "Model Regulation":</b>  | Void            |

**15 Regulatory information**

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Sara**

**Section 355 (extremely hazardous substances):**

None of the ingredient is listed.

**Section 313 (Specific toxic chemical listings):**

None of the ingredients is listed.

- **SARA 311/312:** See Section 2 for the information on physical and health hazards.

**TSCA (Toxic Substances Control Act):**

56-81-5 Glycerol

ACTIVE

25322-68-3 Polyethylene glycol

ACTIVE

**Hazardous Air Pollutants**

None of the ingredients is listed.

(Contd. on page 6)

US

**Safety Data Sheet**  
acc. to OSHA HCS

Printing date 10/27/2023

Version 6

Reviewed on 10/27/2023

**Trade name: PANA VIA 21 ; OXYGUARD II**

(Contd. of page 5)

· **Proposition 65**· **Chemicals known to cause cancer:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

· **Carcinogenicity categories :**· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value)**

None of the ingredients is listed.

**16 Other information**

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **NFPA ratings (scale 0 - 4)**

Health = 0

Fire = 1

Reactivity = 0

· **HMIS-ratings (scale 0 - 4)**

Health = 0

Fire = 1

Reactivity = 0

· **Abbreviations and acronyms:**

DOT: US Department of Transportation

ELINCS: European List of Notified Chemical Substances

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

· **\* Data compared to the previous version altered.**

US