

Safety Data Sheet
according to HPR, Schedule 1

Printing date 01/25/2024



Version 4

Reviewed on 01/25/2024

1 Identification

- **Product identifier**
- **Trade name:** PANAVIA F 2.0 ; PASTE A
- **Article number:** CA089-PA
- **Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the mixture** Dental materials
- **Details of the supplier of the safety data sheet**
- **Supplier**
Kuraray America, Inc.
32 Old Slip, 7th Floor, New York,
NY 10005 U.S.A.
Tel: 800-879-1676
Website: www.kuraraydental.com
- **Further information obtainable from:** Supplier
- **Emergency telephone number:**
CHEMTREC; day and night, within USA and Canada
1-800-424-9300 CCN706984 or +1 703-527-3887 (collect calls accepted)

2 Hazard identification

- **Classification of the substance or mixture**
Skin Sensitizer - Category 1 H317 May cause an allergic skin reaction.
Carcinogenicity - Category 1A H350 May cause cancer.
- **Label elements**
- **GHS label elements**
The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms**
 
GHS07 GHS08
- **Signal word** Danger
- **Hazard-determining components of labeling:**
Quartz (SiO₂)
benzoyl peroxide
- **Hazard statements**
H317 May cause an allergic skin reaction.
H350 May cause cancer.
- **Precautionary statements**
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 If on skin: Wash with plenty of water.
P308+P313 IF exposed or concerned: Get medical advice/attention.
P321 Specific treatment (see on this label).
P362+P364 Take off contaminated clothing and wash it before reuse.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Safety Data Sheet
according to HPR, Schedule 1

Printing date 01/25/2024

Version 4

Reviewed on 01/25/2024

Trade name: PANAVIA F 2.0 ; PASTE A

(Contd. of page 1)

3 Composition/Information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of substances listed below with nonhazardous additions.

· **Dangerous components:**

14808-60-7	Quartz (SiO ₂) ⚠ Carcinogenicity - Category 1A, H350	> 70%
94-36-0	benzoyl peroxide ⚠ Organic Peroxides – Type B, H241; ⚠ Eye Irritation - Category 2A, H319; Skin Sensitizer - Category 1, H317	< 1%

- **Other ingredients:**
Hydrophobic aromatic dimethacrylate
Hydrophobic aliphatic dimethacrylate
10-Methacryloyloxydecyl dihydrogen phosphate
Silanated amorphous silica
dl-Camphorquinone
Initiator

4 First-aid measures

- **Description of first aid measures**
- **General information:**
Avoid contact with eyes.
Avoid contact with skin.
- **After inhalation:**
Supply fresh air and to be sure call for a doctor.
In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If symptoms persist consult doctor.
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Use fire fighting measures that suit the environment.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**
Use personal protection recommended in Section 8. Isolate the hazard area and deny entry to unnecessary and unprotected personnel.
- **Environmental precautions:** No special measures required.
- **Methods and material for containment and cleaning up:**
Dispose contaminated material as waste according to section 13.

(Contd. on page 3)

Safety Data Sheet
according to HPR, Schedule 1

Printing date 01/25/2024

Version 4

Reviewed on 01/25/2024

Trade name: PANAVIA F 2.0 ; PASTE A

(Contd. of page 2)

- Ensure adequate ventilation.
- Contain spill, then place in a suitable container. Minimize dust generation. Do not flush to sewer or allow to enter waterways.
- Vacuum or sweep material and place in a disposal container. Provide ventilation.
- Use appropriate Personal Protective Equipment (PPE).
- **Reference to other sections**
- See Section 7 for information on safe handling.
- See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.

7 Handling and storage

- **Precautions for safe handling**
- Ensure good ventilation/exhaustion at the workplace.
- Open and handle receptacle with care.
- Avoid contact with skin and eyes, Do not swallow. Avoid generating dust. Do not breathe dust/fume/gas/mist/vapors /spray.
- **Information about protection against explosions and fires:** Keep respiratory protective device available.
- **Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store in a refrigerator (2-8°C).
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**
- Keep receptacle tightly sealed.
- Protect from heat and direct sunlight.
- Store in cool, dry conditions in well sealed receptacles.
- **Specific end use(s)** No further relevant information available.

8 Exposure controls/ Personal protection

· **Control parameters**

· **Components with limit values that require monitoring at the workplace:**

14808-60-7 Quartz (SiO₂)

EL	TWA: 0.025 mg/m ³ ACGIH A2; IARC 1
EV	TWA: 0.10* mg/m ³ *respirable fraction

94-36-0 benzoyl peroxide

EL	TWA: 5 mg/m ³
EV	TWA: 5 mg/m ³

- **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
- **Additional information about design of technical systems:** No further data; see section 7.
- **Personal protective equipment:**
- **General protective and hygienic measures:**
- Keep away from foodstuffs, beverages and feed.
- Immediately remove all soiled and contaminated clothing.
- Wash hands before breaks and at the end of work.
- Store protective clothing separately.
- **Breathing equipment:**
- Suitable respiratory protective device recommended.

(Contd. on page 4)

Safety Data Sheet
according to HPR, Schedule 1

Printing date 01/25/2024

Version 4

Reviewed on 01/25/2024

Trade name: PANAVIA F 2.0 ; PASTE A

(Contd. of page 3)

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

· **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact break trough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:** Safety glasses

9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· Physical state	Solid
· Color:	Colored
· Odor:	Odorless
· Melting point/Melting range:	Not determined.
· Boiling point/Boiling range:	Not determined.
· Flammability (solid, gaseous):	Not determined.
· Explosion limits:	
· Lower:	Not determined.
· Upper:	Not determined.
· Flash point:	Not applicable.
· Auto ignition temperature:	Not determined.
· Decomposition temperature:	Not determined.
· pH-value:	Not determined.
· Viscosity:	
· Kinematic viscosity	Not determined.
· Dynamic:	Not determined.
· Solubility in / Miscibility with	
· Water:	Insoluble.
· Partition coefficient (n-octanol/water):	Not determined.
· Vapor pressure:	Not applicable.
· Density at 20 °C:	2.2 g/cm ³
· Relative density	Not determined.
· Relative vapor density	Not determined.

· **Other information**

· **Important information on protection of health and environment, and on safety.**

· **Danger of explosion:** Product does not present an explosion hazard.

(Contd. on page 5)

Safety Data Sheet
according to HPR, Schedule 1

Printing date 01/25/2024

Version 4

Reviewed on 01/25/2024

Trade name: PANAVIA F 2.0 ; PASTE A

(Contd. of page 4)

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
No decomposition if used and stored according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** Carbon monoxide and carbon dioxide

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:** Based on available data, the classification criteria are not met.
- **on the skin:** Based on available data, the classification criteria are not met.
- **on the eye:** Based on available data, the classification criteria are not met.
- **Sensitization:**
May cause an allergic skin reaction.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** May cause cancer.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **Specific target organ toxicity - single exposure**
Based on available data, the classification criteria are not met.
- **Specific target organ toxicity - repeated exposure**
Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **Subacute to chronic toxicity:**

· **Carcinogenic categories**

· IARC (International Agency for Research on Cancer)		
14808-60-7	Quartz (SiO ₂)	1
· NTP (National Toxicology Program)		
14808-60-7	Quartz (SiO ₂)	K

12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects**
- **Additional ecological information:**
- **General notes:** Not hazardous for water.

13 Disposal considerations

- **Recommendation:**
Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

(Contd. on page 6)

Safety Data Sheet
according to HPR, Schedule 1

Printing date 01/25/2024

Version 4

Reviewed on 01/25/2024

Trade name: PANAVIA F 2.0 ; PASTE A

(Contd. of page 5)

- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to official regulations.

14 Transport information

· UN-Number	Void
· DOT/TDG, IMDG, IATA	Void
· UN proper shipping name	Void
· DOT/TDG, IMDG, IATA	Void
· Transport hazard class(es)	Void
· DOT/TDG, IMDG, IATA	Void
· Class	Void
· Packing group	Void
· DOT/TDG, IMDG, IATA	Void
· Environmental hazards:	Not applicable.
· Special precautions for user	Not applicable.
· Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code	Not applicable.
· UN "Model Regulation":	Void

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

- **Canadian substance listings:**

· Canadian Domestic Substances List (DSL)	
14808-60-7	Quartz (SiO ₂)
94-36-0	benzoyl peroxide

· Canadian Non-Domestic Substances List (NDSL)
None of the ingredients is listed.

· Canadian Ingredient Disclosure list (limit 0.1%)
None of the ingredients is listed.

· Canadian Ingredient Disclosure list (limit 1%)
14808-60-7 Quartz (SiO ₂)

- **National regulations:**

- **Information about limitation of use:**

Workers are not allowed to be exposed to the hazardous carcinogenic materials contained in this preparation. Exceptions can be made by the authorities in certain cases.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases**

H241 Heating may cause a fire or explosion.
H317 May cause an allergic skin reaction.

(Contd. on page 7)

Printing date 01/25/2024

Version 4

Reviewed on 01/25/2024

Trade name: PANAVIA F 2.0 ; PASTE A

H319 Causes serious eye irritation.

H350 May cause cancer.

· Abbreviations and acronyms:

DOT: US Department of Transportation

ELINCS: European List of Notified Chemical Substances

· * Data compared to the previous version altered.

(Contd. of page 6)

CA